

varian 

chromatography systems

2700 mitchell drive/walnut creek/ca 94598-1675/u.s.a.



8200/SPME AUTOSAMPLER OPERATOR'S MANUAL

3800 GC Users




Table of Contents

SECTION 1

INTRODUCTION	5
8200 AutoSampler Features.....	5
Injection Methods.....	6
Liquid Injection	6
SPME Injection.....	6
Large Volume Injection.....	6
Ambient Headspace Injection	6
System Description.....	7
Storage Module.....	7
Pneumatics Module.....	9
Injection Tower Module	11
8200 AutoSampler Specifications.....	15
Performance.....	15
Physical	15
Environmental	16

SECTION 2

INSTALLATION	17
Unpack and Inspect the 8200 AutoSampler	18
Unpack the 8200 AutoSampler	18
Inspect the 8200 AutoSampler.....	18
Remove the Carriage Tie Down.....	19
Install the Mounting Plates	19
Connecting the 8200 Cables.....	20
Connecting the 8200 Communications Cable to the 3800 GC	20
Mount the 8200 AutoSampler	21
Connect Power to the 8200 AutoSampler	23
Vial Retainer Alignment	25
Installing the 8200 Star Workstation Driver	25
How to Configure the AutoSampler	26
Specify the 8200 Hardware Configuration.....	27
Final Installation Check.....	29

SECTION 3

LIQUID INJECTION SETUP	31
Connect the Pneumatics.....	31
Install the Solvent Reservoir Bracket.....	32
Route the Solvent and Waste Lines	33
Install the Syringe/Needle Assembly	34
Align the Syringe Needle with the Injector Port.....	36

Build the 8200 Liquid Injection Method on the Star Chromatography Workstation.....38
 Specify the Configuration Window Options39
 Save Your Method45
 Final Liquid Injection Method Checklist.....46

SECTION 4

LIQUID INJECTION OPERATION.....47
 Sampling47
 Loading the Carrousel.....47
 Solvents and Solvent Reservoirs47
 Sample Vials50

SECTION 5

PERFORMING INJECTIONS51
 Injecting a Single Sample51
 Activate a Method51
 Inject the Sample54
 Creating and Running a SampleList.....55
 Build the SampleList56
 Fill in the SampleList.....57
 How to Set the Default Settings58
 How to Delete SampleList Entries59
 Downloading a Method from the Sample List.....60
 Begin the SampleList.....62
 Building and Running a Sequence62
 Open the Active Sequence Window63
 How to Add a Line.....64
 Set the Action, Method and SampleList/Log.....64
 Begin Your Sequence64
 Automix.....65

SECTION 6

SPME INJECTION SETUP69
 Build a New SPME Injection Method.....69
 Specify the Configuration Window Options70
 Save Your Method72
 Install the Fiber Holder on the Syringe Carriage72
 Align the Fiber with the Injector Port74
 SPME Fiber Calibration Procedure.....76
 Adjust the Vial Retainer Arm Height.....77
 Final SPME Injection Method Setup Checklist.....78
 SPME Adjustments79
 SPME Parameters81

SECTION 7

SPME INJECTION OPERATION	83
Sampling	83
Loading the Carrousel	83
Sample Vials	83
Conditioning the Fiber	84
Sample volumes.....	85
Setting the SPME/GC Parameters	85
Running the Test Sample.....	85

SECTION 8

SPME III WITH AGITATION.....	89
Introduction	89
Operation of the SPME Agitation Accessory	90
Principles of Operation	90
Description of the Agitation Device	90
Installation.....	91
Setup and Operation.....	92
Installation of Agitator Septum	93
Installation of SPME Fiber Holder	93
Building a Method that Includes Agitation	94
Changing From SPME III to Liquid Sampling.....	94
Troubleshooting and Maintenance	97
Method Development Hints	97

SECTION 9

AMBIENT HEADSPACE INJECTION SETUP.....	99
Build the Ambient Headspace Injection Method	99
Specify the Configuration Window Options	101
Save Your Method.....	102
Connect the Pneumatics	103
Remove Solvent from the Solvent Lines.....	104
Flush the Solvent Lines	105
Install the 100- μ L Syringe Assembly.....	105
Align the Syringe Needle with the Injector Port.....	106
Final Ambient Headspace Setup Checklist.....	107

SECTION 10

AMBIENT HEADSPACE INJECTION OPERATION	109
Sampling.....	109
Loading the Carrousel	109
100- μ L Syringe Headspace Injection Method	109

SECTION 11
AUTOMATED LARGE VOLUME INJECTION SETUP AND OPERATION111
Introduction.....111
Components Necessary for LVI111
 Conditions for LVI with Various Solvents112

SECTION 12
MAINTENANCE/TROUBLESHOOTING117
Syringe Carriage Guide Rods.....118
Syringe Carriage Drive Belt118
Solvent Compatibility119
Syringe and Syringe Needle Maintenance120
 Syringe Care120
 Needle Care and Replacement120
Replacing the Main Power Fuse121
Instrument Troubleshooting121

SECTION 13
8200 PARTS AND SUPPLIES129

SECTION 14
APPLICATION NOTES.....131

INDEX139

Section 1

Introduction

The Varian 8200 AutoSampler is a low loss, low volume autosampler. This manual refers specifically to the mounting and operation of the 8200 AutoSampler with the 3800 Gas Chromatograph. The AutoSampler is designed to be user installed.

8200 AUTOSAMPLER FEATURES

The 8200 AutoSampler.

- Supports multiple sampling injection methods, including Liquid, Solid Phase MicroExtraction (SPME), and Ambient Headspace Injection.
- Accepts a 10 μ L or 100 μ L syringe.
- Uses a solvent flush sample injection technique (Standard Liquid Injection).
- Injects neat, viscous, and volatile samples (Liquid Injection).
- Accommodates 100, 120, 220, and 240V power.
- Is compatible with the 1079, 1041 and 1061 liquid injectors on the 3800.
- Must be controlled from the Varian Star Chromatography Workstation..
- Permits syringe washing with up to two solvents (Liquid Injection).
- Allows you to program the syringe needle depth in the vial (Liquid or Ambient Headspace or SPME Injection).
- Can be used with a 48-vial (2 mL) or 12-vial (10 mL) carousel.
- Consists of a pneumatics module, injection tower, and storage module in a single compact unit.
- Permits easy access to the pneumatics.
- Supports agitation of the sample during SPME.

INJECTION METHODS

You may perform liquid, Solid Phase MicroExtraction (SPME), or ambient headspace injections with the 8200 AutoSampler.

LIQUID INJECTION

Four liquid injection modes (Standard Solvent Flush, Volatile, Neat, and Viscous) are pre-programmed with optimized parameters for specific samples. A fifth mode, the user-defined mode, allows you to optimize parameters for specific samples. The sandwich, or solvent-flush, injection technique reduces mass discrimination and provides optimum accuracy and precision. It is used in the standard, viscous and volatile modes and can be programmed in the user defined mode.

SPME INJECTION

Solid Phase MicroExtraction (SPME) is an extraction technique for organic compounds in aqueous samples in which analytes are absorbed directly from the sample or the sample headspace onto a fused silica fiber that is coated with an appropriate stationary phase. While the fiber is inserted in the sample or its headspace, the analytes partition from the sample matrix into the stationary phase until equilibrium is reached. The fiber is then inserted into a capillary injector port of the gas chromatograph where it is heated, and the analytes are rapidly thermally desorbed into the capillary GC column for analysis. In the 3800, SPME is typically used with the 1079 injector.

LARGE VOLUME INJECTION

The 8200 can accommodate a 100 μL syringe for automated large volume injection. Injection volumes can be programmed from 1-100 μL .

AMBIENT HEADSPACE INJECTION

Ambient headspace analysis can be used to determine organic solvents in aqueous matrices (e.g., ethanol in blood or urine, solvents in waste water). Generally, volatiles in a solid matrix are not detected in the headspace of an unheated system. A 100- μL syringe is used with the 8200 AutoSampler for ambient headspace analysis. Note that this is not the same syringe as that used for large volume injections.


SYSTEM DESCRIPTION

The 8200 AutoSampler mounts over the injectors on the top of the Varian 3800 Gas Chromatograph. Mounting brackets are included in the 8200 accessory kit which attach to the 3800 injector cover.

A power entry module located on the rear panel of the pneumatics module supplies power to the AutoSampler. The rear panel also contains the main power switch, fuse(s) and voltage selection components.

You can operate the AutoSampler at 100, 120, 220, and 240 Volts, 50/60 Hz, with single phase, or split phase (phase/phase) power. For information about selecting the appropriate voltage for the AutoSampler, see **Connect Power to the 8200 AutoSampler** on Page 23. For information about changing fuses, see **Replacing the Main Power Fuse** on Page 121.

The AutoSampler has three modules: the storage module, pneumatics module, and injection tower module, as shown in the following diagram.



STORAGE MODULE


The storage module consists of a 48- or 12-vial carousel, turntable, position sensor board, vial sensors, and two synchronous motors that rotate the turntable and move it in and out of the sampling position.

The power ON indicator (LED) is located on the front right side of the storage module. This LED lights up when the power cord is properly connected and the AutoSampler power switch is ON. The power switch is located on the rear of the tower.

Carrousels

- 48-vial Carrousel

The 48-vial carrousel is constructed from a chemically resistant polymer. It is 5-3/4 inches in diameter by 1-1/4 inch deep, and can hold up to forty-eight 2-mL vials. A center knob allows you to lift and reposition the turntable. A small raised indicator marks the Home vial position. See the following diagram.




The carrousel label designates the location of the 48 vials, as well as the location of the Home position (Vial 1). You position the label on the carrousel using the keyed cutout. The carrousel label is affixed with a semi-permanent adhesive that permits the label to be easily removed and replaced. Section 13, **8200 Parts and Supplies**, provides ordering information for additional labels.

NOTE: The 48-vial tray can be placed on the turntable in any position.

- 12-vial Carrousel

The 12-vial carrousel is constructed from anodized aluminum. It is 5-3/4 inches diameter by 1-1/4 inches deep, and can hold up to twelve 10-mL vials. A center knob allows you to lift and reposition the turntable. A small raised indicator marks the Home vial position. See the following diagram.




The carousel label designates the location of the 12 vials, as well as the location of the Home position (Vial 1). You position the label on the carousel using the keyed cutout. The carousel label is affixed with a semi-permanent adhesive that permits the label to be easily removed and replaced. Section 13, **8200 Parts and Supplies**, provides ordering information for additional labels.

NOTE: The 12-vial tray must be placed on the turntable with the pin aligned with vial 1.

Turntable

Two synchronous motors drive the turntable: One motor rotates the turntable clockwise, and the other moves it backward or forward. The carousel mounts on the turntable and rests on a vertical pin on the outer edge of the turntable. See diagram below.


NOTE: Make sure the carousel is firmly seated whenever you place it on the turntable.



PNEUMATICS MODULE

The pneumatics module houses the solvent and air solenoid valves that control solvent and air to the solvent reservoirs and pressure regulator; the power transformer; the main power switch, fuses, and voltage selection module; and the Controller PC Board. Gas (at 40 to 60 psig) enters the pressure regulator. The pressure regulator is preset to 20 psi (outlet) at the factory. You should not have to adjust the gas pressure. Gas is distributed to four valves (L1, L2, L3, L5) that control pressure and solvent flow from reservoirs A and B. The solenoid valves in the pneumatics module connect to J8 on the 8200 Controller PC Board.


- A Pressurization Valves
- B Solvent Distribution Valves
- C To Syringe Side Arm
- D Power Entry Module
- E Transformer
- F Gas/Air Input, 40-60 psig
- G Pressure Regulator



Solenoid Valves Assignments

The solenoid valves in the pneumatics module are assigned as follows:

- L1 controls solvent flow from Reservoir B to the syringe side arm.
- L2 controls solvent flow from Reservoir A to the syringe side arm.
- L3 pressurizes solvent Reservoir A.
- L5 pressurizes solvent Reservoir B.



Power Entry Module


The power entry module sits on the lower rear panel of the pneumatics module. It includes a main power switch; main power cord receptacle; fuse(s); and voltage selection card. Before powering up the 8200 AutoSampler, make sure that the white voltage pin is in the hole that matches the voltage rating on the serial number label, as described on Page 23, **Connect Power to the 8200 AutoSampler**. If you need to change a fuse, see Page 121, **Replacing the Main Power Fuse**.

INJECTION TOWER MODULE

The injection tower module houses the syringe carriage and syringe mounting hardware, the syringe carriage drive stepper motor, the waste cup, and the syringe plunger stepping motor drive assembly. The hinged storage module and pneumatics module both bolt to the injector tower module.

The stepping motor and drive belt move the syringe carriage through the vial sampling, injecting, and washing positions. See the following diagram.

- A Syringe Plunger Lead Screw
- B Belt Tensioner/Pulley
- C Syringe Carriage Stop Latch Solenoid
- D To Position #4 on J8
- E Syringe Carriage Drive Belt
- F To Position #7 on J8
- G Syringe Carriage Drive Motor
- H Syringe Carriage




At the end of a run, the syringe carriage returns to the Home position. If power is turned off or if there is a power failure, the stop solenoid opens and the carriage drops to its rest position on the stop latch. When power is restored, the stop solenoid pulls in and the carriage returns to its operating Home position.

Liquid Injection

The liquid injection syringe/needle assembly (P/N 03-918986-00) includes the following:

- 10- μ L syringe with a side inlet fitting
- plunger guide
- non-coring side hole needle (P/N 03-918987-00)


The 10 μL syringe is shown schematically in the following diagram.



- A Needle
- B Syringe Nut
- C Spring
- D Needle Stop
- E PTFE Needle Seal
- F Front Screw Thread
- G Syringe Barrel
- H PTFE Plunger Tip
- I Side Arm Tee
- J Plunger Guide and Thrust Assembly
- K Lock Nut
- L Plunger

SPME Injection

The SPME holder with fiber assembly is shown below.




- A SPME Fiber
- B Septum Piercing Needle
- C Needle Ferrule
- D Retaining Nut
- E Fiber Attachment Needle
- F Barrel
- G Plunger
- H Sealing Septum
- I Color-Coded Screw Hub
- J Slot
- K Retaining Screw

Waste Arm/Cup Assembly

The waste arm, which is used with the Liquid and Ambient Headspace injection methods, is electrically driven to the "out" position by an electrical solenoid; a spring returns it to the "in" position. The waste arm positions the cup under the needle during washing or flushing cycles. The cup remains in the retracted position during sample injection.

When it becomes necessary to change the waste cup septum, replace it only with a Varian septum (P/N 00-997628-00). The Standard 8200 Accessory Kit (03-918704-91) includes a supply of these septa. You install the new septum with the Teflon®-coated surface down, as seen in the diagram below.

There is no fixed schedule for changing the waste cup seal. You should, however, change it if you see signs of leakage or sample carryover.



Vial Retainer Arm


The vial retainer arm is a metal bar clamped in a plastic holder. It is positioned just above the hole in the sample vial cap. The retainer prevents the vial from being lifted out of the carousel when the needle is withdrawn. During injection, as the syringe carriage descends, it moves the retainer out of the syringe carriage path. When the syringe carriage moves up from the inject position, the retainer returns to its normal position over the vial. The following diagram shows the retainer in its correct position.

The vial retainer arm is adjusted at the factory to align with the syringe needle guide. If the vial retainer arm is out of position, the arm could either interfere with the operation of the needle guide, or fail to hold the vial in the carousel during the withdrawal of the syringe needle. Turn the adjusting screw to position the retaining arm as shown.

**CAUTION**

Adjust Height of Vial Retainer Arm When Changing Vial Size.

- A** Vial Retainer Arm Block
- B** Vial Retainer Arm
- C** Syringe Needle Guide Arm
- D** Adjusting Screw



8200 AUTOSAMPLER SPECIFICATIONS

PERFORMANCE

Peak Area: <1.0% RSD typical in standard mode (Liquid Injection) for all injection techniques except isothermal split/splitless (1079) where typical RSD is <1.5%.

Carryover: <0.05% typical with default values.
<0.01% when wash cycle is lengthened by increased viscosity factor or with multiple solvent wash.

PHYSICAL

Height: 19.5 inches (49.5 cm)

Width: 6.0 inches (15.5 cm)

Depth (Including Storage Module): 15.75 inches (40 cm)

Weight: 24 lb (10.9 kg)

Sample Capacity: 48 samples contained in 2-mL screw-cap or crimp-top vials, standard microvials with internally machined conical volume or microvial inserts. These types can be mixed within the sample carousel. Wash vials are not required. Alternately, 12 samples contained in 10-mL sample-volume crimp-top vials.

Liquid Injection Syringe:	Flow-through design to allow dynamic washing; standard 10- μ L syringe with Teflon [®] -tipped plunger. Options include 100- μ L syringe and SPME fiber holder.
Sampling Methods:	Liquid Injection (Standard, Volatile, Neat, Viscous, and User Defined modes); SPME Injection; Large Volume Injection and Ambient Headspace Injection. Note that when controlling the 8200 locally from the 3800, the only modes available are the Standard, Volatile, Neat and Viscous modes.
Programmable Uptake Speed	1.0 to 5 μ L/sec
Programmable Injection Speed:	0.2 to 10 μ L/sec
Programmable Residence Time:	0 to 10 minutes
Programmable Hot Needle Time:	0 to 1.0 minute
Sample Size:	Programmable from 0.1 to 10.0 μ L with 10- μ L syringe. Programmable from 1 μ L to 100.0 μ L with 100- μ L syringe.
Sample Waste:	Less than 0.1 μ L per injection in standard solvent-flush mode (Liquid Injection).
Minimal Sample Volume (Liquid Sampling):	2-mL vials: 300 μ L Microvials: 10 μ L 10-mL vials: 5 mL
External Wash Solvent:	Choice of solvent A and/or solvent B for Liquid Injection.
Sample Coding:	Vials 1-48, or vials 1-12.

ENVIRONMENTAL

For optimum 8200 AutoSampler performance, check that the operating environment meets the following specifications:

Temperature:	Operating: 10 to 40°C Storage: -20 to 65°C
Humidity:	Operating: to 80% non-condensing
Power:	90 to 264 Vac operating range with four selectable voltages: 100, 120, 220, 240 Vac (\pm 10%) Frequency: 50/60 Hz Power Consumption (Maximum): 55 VA

Section 2

Installation

To install the 8200 AutoSampler prior to setting up an injection method, you need to:

- Unpack and inspect the 8200 AutoSampler.
- Install the mounting hardware on the 3800 GC injector cover.
- Install the SPME hardware (SPME Injection only).
- Mount and secure the AutoSampler to the mounting plates.
- Check that the operating voltage label and voltage selection pin agree with the on-site power rating, that a properly rated fuse has been installed, and that the power cord has been connected.
- Install and configure the AutoSampler software.

Except for installing the mounting hardware, which may require up to 15 minutes to complete, none of these tasks should require more than 5 minutes.

NOTE: The 8200 AutoSampler can be programmed from the local 3800 keyboard or from the Star Chromatography Workstation. For local 3800 control of the 8200 AutoSampler, refer to the 3800 Operator's Manual. Note that when controlling the 8200 locally from the 3800, the only modes available are the Standard, Volatile, Neat and Viscous modes. Star Chromatography Workstation communication to the AutoSampler is via the Ethernet® link between the PC and the 3800. Before installing the 8200 ensure the Star Chromatography Workstation is installed and is communicating with the 3800 GC.

Before you open the carton your 8200 AutoSampler arrived in, check for any evidence of damage, e.g., crushed corners, forklift punctures, water stains. If you detect any damage, notify Varian immediately at

Varian Chromatography Systems
2700 Mitchell Drive
Walnut Creek, California 94598-1675
Attention: Manager of Customer Service
Phone (510) 939-2400

or contact your local Varian Sales/Service Center.

UNPACK AND INSPECT THE 8200 AUTOSAMPLER

UNPACK THE 8200 AUTOSAMPLER

Your 8200 AutoSampler will have arrived in a single shipping carton. Arrows indicate the top of the carton.

To unpack the AutoSampler, proceed as follows:

1. Place the carton on the floor near the GC on which you plan to install the AutoSampler. Make sure the carton is top-side up.
2. Cut the two plastic shipping straps and lift the top cover. Remove the inner cover with its attached foam support blocks.
3. Use the two hand holes to lift the inner liner and its attached support blocks from the carton.

NOTE: Save the shipping carton in case you need to ship the instrument back to the factory.

4. Open the inner carton and check that it contains the Standard 8200 Accessory Kit (P/N 03-918704-91). Check the contents of the Standard Accessory Kit against the enclosed list. If there are discrepancies or items missing, notify Varian immediately.
5. Lift the AutoSampler up and out of the shipping carton and place it upright on the bench next to the GC.



CAUTION

Do not lift the AutoSampler by the storage module. Rather, lift at the upper left and lower right corners of the AutoSampler to avoid damaging the storage module hinge.

NOTE: Do not overlook the syringe that accompanies your AutoSampler. The syringe is packed in a separate small container in the shipping carton beneath the AutoSampler.

INSPECT THE 8200 AUTOSAMPLER

Carefully inspect the 8200 AutoSampler for damage. If you find any damage, or if you find that parts are missing, notify Varian immediately.

REMOVE THE CARRIAGE TIE DOWN

To remove the carriage tie down, proceed as follows:

1. While holding the syringe carriage, cut the plastic tie wrap that holds the syringe carriage to the top cover.
2. Lower the syringe carriage to the stop.
3. Carefully remove the tie wrap from the AutoSampler.

INSTALL THE MOUNTING PLATES


To install the mounting plate(s), proceed as follows:

1. Set the mounting plate(s) on the injector cover so that the screw holes align with the injector of choice.
2. Secure the plate(s) to the injector cover with two 8-32 screws and washers. The screws go through the slightly elongated openings in the plate. The elongation permits you to position the AutoSampler over the injector nut.

Although mounting plates can be moved, ideally one plate should be installed per injector.

8200 Injection Module

A Short Plate



CONNECTING THE 8200 CABLES

Prior to installing the 8200 AutoSampler, you must connect its communication cable to the 3800 GC.

CONNECTING THE 8200 COMMUNICATIONS CABLE TO THE 3800 GC


- To install this hardware, you will need a **Phillips** screwdriver:



**WARNING:
SHOCK HAZARD**

**Dangerous Voltages Exposed When High Voltage Cover is removed.
Unplug Power Cord.**

1. Turn off the power to both the 3800 GC and the AutoSampler. The Varian 3800 GC power switch is located on the top rear of the GC. The AutoSampler switch is located at the lower left rear of the AutoSampler.
2. The 8200 communications cable is already installed in the 8200. Locate this cable identified by a 25-pin connector on the end.
3. Remove the top left cover from the 3800 GC and locate the (J4) SID1 25-pin connector on the exposed top panel of the GC.
4. Run the 8200 cable through the oblong slot on the rear of the GC and connect the cable to the (J4) SID1 connector.



Be sure to use the mounted plastic clips to keep air and solvent lines away from hot vent.



WARNING

Melting of air and solvent lines may result. Route air and solvent lines away from hot vent using mounted plastic clips.


5. Replace the top left cover on the GC.

MOUNT THE 8200 AUTOSAMPLER

To mount the 8200 AutoSampler on the GC, proceed as follows:

1. Position the AutoSampler over the selected injector nut.

2. Carefully slide the AutoSampler on the mounting plate until the ball plunger screws drop into the ball plunger detents in the mounting plate.
3. With the storage module swung open for access, press down and engage the retractable captive screws in the injection tower base.
4. Turn the captive screws clockwise to secure the AutoSampler in place. See below.



The spring-loaded balls in the ball plunger screws should rest in detents in the mounting plates. They serve to accurately position the AutoSampler. The ball plunger screws were preset in the chassis at the factory and should not require further adjustment. However, you can readjust them if the balls do not engage the detent properly, or if the AutoSampler rocks on its base after you have tightened the captive screws.


If you need to adjust the ball plunger screws, use a **thin, flat blade screwdriver**.

To adjust the ball plunger screws, proceed as follows:

1. Turn the captive screws slightly (clockwise or counterclockwise) with the thin flat blade screwdriver.
2. Test the adjustment by loosening the captive screws and sliding the AutoSampler in and out of the ball plunger detents. When properly adjusted, the AutoSampler will move smoothly on the mounting plate. Four threaded nylon feet in the injection tower base provide easy sliding for positioning, and prevent scratching of painted or finished surfaces.
3. The ball plunger screws extend 0.055 inches from bottom of the base of the injection tower.
4. Retighten the captive screws.


7. Orient the selector card so that you can read the voltage to be used at the bottom (the indicator pin points out).
8. Insert the voltage selector card into the housing with the desired voltage showing on the edge. The printed side of the selector card should face upwards toward the IEC connector.
9. Replace the cover and check that the pin is showing in the hole of the desired operating voltage.
10. Connect the power cord to J18 on the panel of the AutoSampler. Plug the other end of the power cord into a wall receptacle.

NOTE: If you intend to operate the AutoSampler at a specific voltage for an extended period of time, write the operating voltage on the label on the rear panel of the AutoSampler.



VIAL RETAINER ALIGNMENT

- A Vial Retainer Arm Block
- B Vial Retainer Arm
- C Syringe Needle Guide Arm
- D Adjusting Screw



Follow Steps 1 through 5 on Page 36, **Align the Syringe Needle with the Injector Port.**

1. Lower the carriage by hand until the syringe needle guide arm is level with the vial retainer arm
2. Adjust setscrew so that the retainer arm just clears the guide arm.
3. Raise carriage to upper limit of travel, then lower until it rests on the stop solenoid.
4. Close the storage module.

INSTALLING THE 8200 STAR WORKSTATION DRIVER

To install the 8200 AutoSampler software as an add-on to the Star Chromatography Workstation, proceed as follows:

1. Verify that version 4.51 (or later) of the Star Chromatography Software is installed on your computer. You can do this by selecting **About...** under **Help** in the menu bar of any of the Star Workstation applications and checking the version number. If you are running an earlier version, install version 4.51 (or later) before proceeding.
2. Close all Star Workstation Applications, then from the Start button click on the Module Installation icon in the Varian Star 4.5 folder. This brings up a dialog box showing a list of the currently installed drivers for the Star Workstation.
3. To add or update the 8200 driver, press the **Add...** button. This displays the Add Module Drivers dialog box.

4. Click on Floppy Drive/Directory. Specify the location from which you want to install the driver. By default this location is A:\. If the floppy drive used for the installation is not A:, substitute the appropriate designation (for example, B:\).
5. Insert the floppy disk for 8200/SPME Driver Software (P/N 03-910729-00) in the specified floppy drive, then press the **OK** button. The Add Module Drivers dialog box will now list the drivers available for installation.
6. Click on **Install All**. An asterisk will be added in front of the driver to indicate that it has been installed.
7. Click on **Done**. You are returned to the Varian Star Workstation Drivers window. Verify that the 8200 driver has been added to the list of installed drivers in the dialog box.
8. Click on **Done**. This completes driver installation. You can use the Module Installation program any time to add or remove drivers from your Workstation. For more details on the program's operation, see the Workstation's On-Line Help files.

HOW TO CONFIGURE THE AUTOSAMPLER

Once you have installed the 8200 AutoSampler software, the first thing you will want to do is to configure your GC for AutoSampler operation. You do this from the AutoSampler's System Control Configuration page.

NOTE: Before you configure the software for AutoSampler operation, be sure you have turned on the AutoSampler. The power switch is located on the rear of the tower.

To configure the AutoSampler, proceed as follows:

1. From the Start page click on the System Control/Automations button or from the Windows Start button, then click on the System Control icon in the Varian Star 4.5 folder.

The System Control Configuration window appears. The Configuration window accommodates up to four AutoSamplers, each connected to a different GC, at any given time. A quadrant of the System Control Configuration screen is reserved for each of the four instruments, i.e., Instrument 1, Instrument 2, Instrument 3, and Instrument 4.


NOTE: Until you have configured an Instrument, a Help screen and Setup COMM Ports screen will appear each time System Control is started. Close these screens to continue.

When each 8200 connects, it will be represented by an icon at the bottom of the Configuration Window in the Available Modules area.

The icon will have a label reading nn:xxx where nn is the address that corresponds to the 3800 GC to which the 8200 is connected and xxx is the host name of the corresponding 3800 GC. The addresses are assigned as follows:


3800 Address	44	45	46	47
8200 Bus Address	26	30	34	38

1. Click and drag the 8200 AutoSampler icon from Available Modules into the AutoControl Module of the Instrument workspace where its corresponding 3800 GC is configured.




SPECIFY THE 8200 HARDWARE CONFIGURATION

1. Double-click in the quadrant of the System Control Configuration window corresponding to the instrument whose hardware configuration you wish to specify. Double-click on icon in the lower left corner. The following window appears.



2. Click on the Hardware button as soon as it becomes active. The following window appears.



3. Under Sampling Device, select the sampling technique you wish to use, i.e., 10- μ L syringe, 100- μ L syringe, or SPME device.
4. There are two check boxes that configure the synchronization signals for different combinations of hardware. These boxes must be configured as follows:

AutoSampler Mode	Instrument	“Wait for Ready In”	“Invert Ready In”
SPME or Liquid Sampling	3800 <i>with</i> Star Workstation Control	Unchecked	Unchecked
SPME or Liquid Sampling	3800 <i>with</i> Saturn MS	Checked	Checked

FINAL INSTALLATION CHECK

Installation of the 8200 AutoSampler is now complete. Before setting up an injection method for the AutoSampler, make sure you have

- Installed the mounting plate and aligned and secured the autosampler cover for the desired injector.
- Installed the communications cable to the 3800 GC.
- Checked that the operating voltage label and voltage selection pin agree with the on-site power rating, that a properly rated fuse has been installed, and that the power cord has been connected.

Section 3

Liquid Injection Setup

To set up the 8200 AutoSampler for Liquid Injection, you will

- Connect the pneumatics.
- Install the solvent reservoir bracket.
- Route the solvent and waste lines.
- Install the syringe/needle assembly.
- Align the syringe needle with the injector port.
- Build a liquid injection method.

With the possible exception of connecting the pneumatics, none of these individual tasks should take you more than five minutes to complete.

NOTE: **Liquid Sampling requires removal of the SPME agitation accessory. If this device is installed on your autosampler refer to Changing From SPME III to Liquid Sampling on Page 94.**

CONNECT THE PNEUMATICS

The 8200 AutoSampler requires clean, chromatographically pure compressed (40 to 60 psig) air or nitrogen to operate. Because the gas contacts the solvent, a chromatographic grade gas is recommended. It is not necessary to use a separate air or nitrogen tank for the AutoSampler; nitrogen or air tanks that are currently used for detectors may have a 1/8-inch Tee fitting (P/N 28-691944-00) installed for pressurizing the AutoSampler.

NOTE: **If you have already connected the pneumatics, proceed to Install the Solvent Reservoir Bracket on Page 32.**


To connect the pneumatics, proceed as follows:




**WARNING:
EYE HAZARD**

Wear Proper Eye Protection. To prevent possible injury from a flailing gas line, shut off the gas supply valve before disconnecting the line.

1. Connect the supply cylinder to the “barbed” fitting on the rear panel of the AutoSampler using the 1/8-inch OD polyurethane tubing from the Standard Accessory Kit. The barbed fitting is labeled GAS INLET 40-60 psig.



2. To make the connection at the GAS INLET, press the 1/8-inch OD tubing onto the barbed hose connector. Tug firmly on the tubing to make sure it is fully gripped on the fitting.



Detail of Gas Quick Connect Fitting


INSTALL THE SOLVENT RESERVOIR BRACKET

The solvent reservoir bracket holds two 180-mL solvent bottles and is secured to the left side panel of the 3800 GC.

To install the solvent reservoir bracket, you will need the following:

- 2.5-inch strip of adhesive-backed neoprene tape. The tape comes with the Standard Accessory Kit (P/N 03-918704-91).
 - No. 2 Phillips screwdriver.
1. Remove the upper left corner screw on the 3800 side panel.
 2. Affix the 2.5-inch strip of neoprene tape to the rear surface of the reservoir bracket. Do not cover the bracket center hole.

3. Align the center hole of the reservoir bracket with the side panel screw hole. Secure the bracket to the GC with the panel screw.



ROUTE THE SOLVENT AND WASTE LINES

Two sets of reservoir cap/line assemblies connect to their assigned solenoid valves in the pneumatics module. The reservoir cap/line assemblies were wrapped in plastic bags at the factory.

To route the solvent and waste lines, proceed as follows:

1. Before installing the reservoirs, unwrap the assemblies and check them for loose fittings.
2. The reservoir cap/line assemblies are tagged A and B. The 180-mL solvent bottles included in the Standard Accessory Kit are not labeled. To avoid later confusion, clearly label the bottles A and B. Check that the cap assembly A goes to solvent reservoir A and that cap assembly B goes to solvent reservoir B.

Be sure to use the mounted plastic clips to keep air and solvent lines away from hot vent.




WARNING

Melting of air and solvent lines may result. Route air and solvent lines away from hot vent using mounted plastic clips.

3. Route the Teflon® waste line attached to the AutoSampler to a waste bottle.

- To ensure an adequate drain from the solvent reservoir bottles to the waste bottle, position the waste bottle below the level of the waste cup. Avoid introducing kinks or low spots in the drain line.




INSTALL THE SYRINGE/NEEDLE ASSEMBLY

To install the syringe/needle assembly on the syringe carriage, proceed as follows:

- Open the front door of the syringe mounting assembly by simultaneously and firmly lifting *both* side-locking tabs. Raise the tabs to the upper position. This may require considerable force, so hold the carriage assembly down with one hand. Alternatively, you can hold the mount down with both thumbs while you lift the tabs. The locking mechanism will become less stiff with use.

- A Solvent Line
- B Plunger Drive
- C Carriage
- D Side Locking Tabs (2 places)
(Down: door closed; red tab not visible)
- E Syringe Mount
- F Carriage Bracket
- G Needle Guide
- H Vial Retainer Arm
- I Red Tab
- J Syringe Mount Assembly (door open)



NOTE: A red tab is visible above the syringe mount assembly when the door is unlocked.

2. Referring to the diagram, guide the syringe needle into and through the small hole in the needle guide.
 - Press the syringe nut into the slot so that it rests on the upper surface of the syringe carriage bracket.
 - Press the syringe plunger button into the slot of the lead screw nut.
3. Close the syringe mount door.

NOTE: Check that the door is fully closed before pushing the side-locking tabs down. Then, holding the door closed and supporting the carriage with one hand, press the side locking tabs down with the other hand. This locks the syringe into position so that the syringe hex nut does not move. When the tabs are down, the door is locked and the red signal tab is not visible.

4. Connect the solvent line to the syringe side arm. Refer to the diagram.



**WARNING:
EYE HAZARD**

Wear Proper Eye Protection.

The Teflon® tube may disconnect from the syringe and spray solvent.


- Slide the knurled nut back from the end of the line.
 - Press the solvent line firmly against the machined shoulder of the side arm port.
5. Slide the nut up to the port and thread it onto the side port. Tighten the nut finger tight. Do not use pliers.
 6. If you notice a leak during operation of the AutoSampler, further tighten the nut with your fingers.

ALIGN THE SYRINGE NEEDLE WITH THE INJECTOR PORT


NOTE: When you initially mount the AutoSampler on the GC, the alignment of the syringe needle with the needle port of the injector nut may not be precise.

To align the syringe needle with the injector port, proceed as follows:

1. Check that the AutoSampler power switch is OFF. The power switch is located on the rear of the tower.



2. Press down on the storage module latch release (see below) and swing the module out to access the GC injector. When the storage module is swung out (with power off), a contact in the latch assembly activates that turns the syringe carriage stop solenoid on or off each time the carriage position sensor is interrupted.




NOTE: Press the latch release down when closing the storage module.

3. Turn the GC power ON.
4. Turn the AutoSampler power ON.
5. Raise the syringe carriage from its rest position to its top limit. The stop solenoid toggles (you will hear a “click”) to allow you to lower the carriage by hand.



Take care that the carriage does not drop and damage the needle assembly.

6. Loosen the mounting plate screws (see below) so that the AutoSampler can be moved on the GC cover.
 - Lower the syringe carriage by hand until the needle tip is just above the GC injector needle port.



7. Carefully adjust the position of the AutoSampler to bring the syringe needle into precise alignment with the injector needle port.
 - Securely tighten the mounting plate screws
8. Lower the syringe carriage and puncture the injector septum two or three times to ensure that the syringe and syringe needle are properly aligned with the injector port.
9. To return the syringe carriage to its rest position, raise the syringe carriage to its upper limit of travel. Then lower the carriage a fraction of an inch until it engages.
10. Close the storage module.

BUILD THE 8200 LIQUID INJECTION METHOD ON THE STAR CHROMATOGRAPHY WORKSTATION


NOTE: Be sure you have selected the 10- μ L syringe setting in the 8200 Hardware Configuration Window (see Specify the 8200 Hardware Configuration, Page 27).

The method is a record of all the instructions and parameters required to analyze your sample. The 8200 AutoSampler software allows you to create separate methods for each of your sample types.


To build a new Liquid Injection Method, proceed as follows:

1. From the Start page click on the Method Editor button, or from the Windows Start button click on the Method Editor icon in the Varian Star 4.5 folder.

The Method Editor window appears with the Create/Open Method File in the foreground.



2. Select "Create a new method file" and click **OK**.
3. The New Method Section window appears. Select "8200 Stand-Alone AutoSampler" in the Instrument Control section of the window.



- Referring to the following table, use the up and down arrows to set the Bus Address. Click on **OK**.

If the 3800 GC Address is...	Set the Bus Address
44	26
45	30
46	34
47	38

- The New Method Section window closes and the 8200 Configuration icon appears in the lower left hand corner of the Method Editor window.
- Double-click on the 8200 Configuration icon. The Configuration window appears.

SPECIFY THE CONFIGURATION WINDOW OPTIONS

- Click on the downward-facing arrow to select the Carrousel Type: 48 Vials or 12 Vials.

The Configuration window is titled "Configuration" and contains the following sections:

- Carrousel Type:** A dropdown menu currently set to "48 Vials".
- AutoSampler Timing:**
 - EndTime: A numeric input field set to "1.00" with up and down arrows.
 - Prep Ahead
- Sampling Modes:**
 - Liquid:**
 - Standard Mode
 - Volatile Sample
 - Neat Sample
 - Viscous Sample
 - User Defined
 - SPME Mode
- Liquid Sampling Options:**
 - Solvent A Wash: A dropdown menu.
 - 10 uL Syringe: A dropdown menu.

At the bottom of the window are "Save" and "Cancel" buttons.

- If Prep Ahead is checked, the AutoSampler will start preparing for the next injection before the GC is ready. This feature can be used to minimize the wasted time between runs. If this feature is used, however, the GC End Time should be entered. If a GC cycle time is entered that is shorter than the actual time, the sample will be held in the syringe before injection; this may result in loss of precision.

If you choose not to select Prep Ahead, leave the End Time setting at the default setting of 1.00 minute or enter the value directly.

3. Under AutoSampler Timing, set the EndTime (minutes) using the up and down arrows to change the default value of 1.00 minute.
4. Select the sampling mode type under Liquid in the Sampling Modes section of the Configuration window. There are five sampling modes to choose from. The Standard, Volatile, Neat, and Viscous modes are pre-programmed with optimized parameters for specific samples. The fifth liquid mode, User-Defined, allows you to optimize the parameters yourself.

- **Standard Sample Mode.** The standard sampling mode uses the solvent flush injection technique (Page 47, **Solvents and Solvent Reservoirs**). Good results are almost always obtained using this mode, so you should try it first if you are developing a new method.

The standard sampling mode is a solvent-flush, or sandwich method of sampling and injection. Parameters for this mode of sampling and their default values appear in the 8200 Quick Reference Guide.

- **Volatile Sample Mode.** Select the volatile sample mode when using samples with high vapor pressures at room temperature (e.g., methylene chloride, carbon disulfide, pentane).

To improve the performance of the AutoSampler when it is sampling volatile samples, some sampling parameters are automatically changed to the following:

- AutoSampler Injection Rate: 1.0 $\mu\text{L}/\text{sec}$
 - Upper Air Gap: Off
 - Sample Uptake Speed: 1.0 $\mu\text{L}/\text{sec}$
 - Pause Time: 6 seconds
- **Neat Sample Mode.** Select the neat sample mode when the wash solvent peaks would interfere with the chromatography of the sample. In the neat sampling mode, the sample is drawn into the syringe and then expelled or injected (suck-and-squirt) without the use of air or solvent gaps.

The Neat sample mode is also enabled when User Defined is selected and No is entered for Solvent Flush Sampling. As in the Neat mode, the sample is drawn into the syringe directly without air or solvent gaps and then injected; however when the User Defined mode is selected, you can more carefully tailor the sampling technique to your specific application.

The Neat Sampling sequence proceeds as follows:

- a) A sample (10 μL) is withdrawn from the sample vial at 5 $\mu\text{L}/\text{sec}$. The sample is expelled into the waste arm at 10 $\mu\text{L}/\text{sec}$. This step is repeated once.
 - b) A sample (10 μL) is withdrawn from the sample vial at 5 $\mu\text{L}/\text{sec}$. The sample is expelled (returned) to the sample vial at 10 $\mu\text{L}/\text{sec}$. This step is repeated six times.
 - c) A sample (10 μL) is withdrawn from the sample vial at 1 $\mu\text{L}/\text{sec}$. The sample is expelled (returned) to the sample vial at 10 $\mu\text{L}/\text{sec}$. This step is repeated three times.
 - d) A sample (5 μL) is withdrawn from the sample vial at 1 $\mu\text{L}/\text{sec}$. The lower air gap is formed. An upper air gap is not formed.
 - e) A specified amount of sample (Default = 1.0 μL) is injected into the GC injection port at 5 $\mu\text{L}/\text{sec}$.
 - f) The excess sample is expelled into the waste arm either before or during the syringe washing cycle.
- **Viscous Sample Mode.** Select the viscous sample mode when injecting samples with a higher viscosity than that of water, e.g., ethylene glycol.

Configuration

Carrousel Type: 48 Vials

AutoSampler Timing

EndTime: 1.00

Prep Ahead

Sampling Modes

Liquid

Standard Mode Viscous Sample

Volatile Sample User Defined

Neat Sample

SPME Mode

Liquid Sampling Options

Solvent A Wash

10 uL Syringe

Save Cancel

To improve the performance of the AutoSampler when it is sampling viscous samples, some sampling parameters are automatically changed to the following:

- AutoSampler Injection Rate: 1.0 $\mu\text{L}/\text{sec}$
 - Solvent Wash Time: 40 seconds
 - Needle Residence Time: 0.2 min.
 - Sample Uptake Speed: 1.0 $\mu\text{L}/\text{sec}$
 - Pause Time: 10 seconds
- **User Defined Mode.** Select the user mode when you want to more carefully tailor the sampling technique to your specific application.

The User Defined mode allows you to adjust the following parameters. Parameter defaults and ranges (in brackets) are also identified.

- **Solvent Flush Sampling** YES [No For Neat Sampling Mode]
- **Syringe Wash Time** (sec) 20 [5 - 180]
- **Air Dry After Wash** No [Yes For Air Dry]. May be set only when solvent flush is not used.

Select **Yes** if you want the syringe air dried following the WASH cycle. The duration of the air dry period is fixed by the syringe wash time.


NOTE: Frequent air drying shortens the lifetime of the syringe.

Solvent Plug Size (μL) 1.0 [0.1 to 3.0 or OFF]. This option specifies the size of the solvent plug used for a “solvent flush” sampling technique. The solvent in reservoir A is used to form the solvent plug.

Set the Solvent Plug Size to OFF if you do not want to have a solvent plug drawn into the syringe.

Vial Needle Depth (%) 90 [0-100]. This option determines the depth which the needle tip reaches into the vial. When set to 0, the needle just penetrates the septum. When set to 100, the needle tip almost touches the bottom of the standard 2-mL vial. Needle depths at 0 and 90 were determined with the internally machined glass vial (P/N 66-000121-00, conical internal construction). At 90%, the needle just reaches the bottom of the machined glass vial. Needle position varies slightly for different vial configurations.

Penetration of the needle into the standard 2-mL vial (03-949835-00) for various needle depths is shown below. This diagram also describes needle depths for the machined microvial.



NOTE: Do not select the 100% setting when using microvials. A 100% setting causes the needle to bottom out in the vial and may damage the needle.

Adjust the vial needle depth only if you are using non-standard microvials, or if you must sample from the upper phase of a two-phase system. To compare the depths of non-standard microvials with a sample of the approved internally machined conical microvial, stand the two vials side-by-side on a level surface. If the bottom of the non-standard vial is lower than the bottom of the machined vial, increase the needle depth for very small samples. If you increase this value more than necessary, the needle may hit the bottom of the vial and damage the needle.

When sampling liquid/solid two phase systems, take care that the needle does not become clogged by solid particles. Particulates from the solid phase can plug the needle.

- **Uptake Speed** ($\mu\text{L}/\text{sec}$) 5.0 $\mu\text{L}/\text{sec}$ [1.0 to 5.0]. This option determines the rate at which the plunger is withdrawn while the syringe needle is in the sample vial. The greater the viscosity of the sample, the lower the sample uptake speed used.
- **Upper Air Gap** Yes [No For None]. May only be set if the solvent flush sampling option has been selected.
- This option determines whether or not an upper air gap is included in the solvent flush sample injection. The upper air gap volume is 0.5 μL .
- **Lower Air Gap** Yes [No For None]. May only be set if the solvent flush sampling option has been selected.

This option determines whether or not a lower air gap is included in the solvent flush sample injection. The lower air gap volume is 0.8 μL .

To reduce discrimination among wide boiling point components when using a vaporizing injector, eliminate the upper air gap and increase the size of the solvent plug.


You must include a lower air gap when using the hot needle injection technique. In this injection technique, the empty needle is heated in the injector before the sample is injected. This method provides more efficient vaporization of the sample.

- **Pause Time** (sec) 2 [1-10]. This option introduces a pause between the time the sample is completely drawn into the syringe and the time the syringe needle is withdrawn from the sample vial and injected. A longer pause time improves precision with viscous samples.
- **Hot Needle Time** (min) 0.00 [0.00 to 1.00]. This option specifies the length of time the needle resides in the hot injector before the sample is expelled from the syringe.
- **Injection Rate** ($\mu\text{L}/\text{sec}$) 5.0 [0.2 to 10.0]. The injection rate is the time used to expel the syringe contents into the injector. With sample sizes of 1 μL or less, the default injection rate (5.0 $\mu\text{L}/\text{sec}$) is usually optimum regardless of injector type. With larger injection volumes (>2 mL), use a slower injection rate to prevent flashback in vaporizing injectors. With non-vaporizing injectors, there is an advantage to a slow injection rate when doing solute focusing. To prevent broadening of early eluting peaks with split injection, use a fast injection rate of 10 $\mu\text{L}/\text{sec}$.

- **Needle Residence Time** (min) 0.00 [0.00 to 10.0]. This option specifies the length of time the needle remains in the injector after the sample has been expelled from the syringe. The residence time differs from the hot needle time; the residence time refers to the period that follows the injection of the sample. Include a needle residence time in the sample injection sequence with vaporizing injectors to minimize loss of volatiles.

SAVE YOUR METHOD

1. Once you have specified your Liquid Injection Method settings, click on **Save**. The Configuration window closes.
2. Close the Method Editor window. A window appears that asks you whether you want to save the changes you have made to the method.
3. Click **Yes**. The Save Method File As window appears prompting you to name the new method.



4. Assign a name of up to eight letters and/or numbers to the method, and click **OK**. The .mth extension will be added automatically.

FINAL LIQUID INJECTION METHOD CHECKLIST

Setup of the Liquid Injection Method for your 8200 AutoSampler is now complete. Before operating the AutoSampler as described in Section 4, make sure you have

- Connected the pneumatics.
- Installed the solvent reservoir bracket.
- Routed the solvent and waste lines.
- Installed the syringe/needle assembly.
- Aligned the syringe needle with the injector port.
- Built a liquid injection method

Section 4

Liquid Injection Operation

This section gives additional information that you will need to operate the 8200 AutoSampler using Liquid Injection. After familiarizing yourself with the information in this section, proceed directly to Section 5, **Performing Injections**.

SAMPLING

LOADING THE CARROUSEL

The carousel rotates clockwise and the vial position numbers increase counter-clockwise around the carousel. Observe this order when loading sample vials. Load sample vials according to the vial numbers shown on the carousel label.

The AutoSampler ignores vials without a septum. Seal all vials in the carousel that are to be sampled with a cap and septum.

SOLVENTS AND SOLVENT RESERVOIRS


Your choice of wash solvent and sampling solvent depends on your particular chromatographic application. Certain chromatographic applications may require two solvents of very different volatilities, viscosities or polarities. Or, you may want to select two wash solvents to ensure the sample is adequately washed from the syringe between injections. The objective is to find the best fit between solvent and sample for your particular application.

- **Filling the Solvent Reservoirs.** The solvent reservoir rack attached to the side of the GC holds two 180-mL solvent reservoirs. The sealed reservoir caps hold two lines. One line is the solvent line and leads from the bottom of the reservoir to the syringe. The other shorter line ends at the neck of the reservoir and controls the pressure to the reservoirs. Keep this line above the liquid level in the reservoirs at all times.



**WARNING:
EXPLOSION HAZARD**


Contents Under Pressure. Keep Bottles in Bracket and Liquid Level Below Vent Tube.



NOTE: When changing the solvent in the wash bottle from one solvent to another, use a new, dry filter. The porous filter absorbs solvent which cannot be completely removed by purging. After you change the solvents, purge the system: set the wash time to 180 seconds and repeat the wash cycle as many times as necessary.

- **Leaks in the Solvent Reservoirs.** When properly tightened by hand, the solvent reservoir caps and seals seal adequately. However, if you smell solvent or hear gas escaping from the solvent reservoirs, check to see if the cap is loose or cracked. If necessary, replace the cap and seal.
- **Venting the Solvent Reservoir.** Volatile solvents or toxic solvents must not be vented directly into the atmosphere. Under such conditions, vent the solvent reservoirs into exhaust ducts or fume hoods. Review your laboratory guidelines for handling and venting toxic and highly volatile solvents established by your department or organization.
- **Solvent Flush Mode Sample Injection.** The solvent flush or *sandwich* technique is the standard mode of injection used in the 8200 AutoSampler sampling sequence. This method of syringe loading is used unless you specify otherwise.

Refer to the following diagram for a description of the traditional mode versus the solvent flush mode of syringe loading. Once you have entered the values for the solvent plug and sample size, this mode of syringe loading is continued automatically until changed.



The solvent flush syringe loading sequence proceeds as follows:

1. The waste arm moves to the out position (over the injector). The syringe carriage moves down to insert the needle into the waste cup.
2. The syringe plunger withdraws from the syringe to its upper limit. Solvent flows from the side port into the syringe barrel and to the waste cup. Solvent flows through the syringe for a length of time equal to one-half the wash time.
3. The plunger moves down to the value entered for the solvent plug size. If you entered 2 μL , the plunger moves to the 2 μL mark on the syringe barrel.
4. The syringe carriage moves to its top position. In some methods, the syringe plunger moves up to introduce the upper air gap.
5. The carousel moves a sample vial under the needle. The syringe carriage moves down and the syringe needle pierces the septum and enters the sample vial. The syringe plunger moves up to draw a specified sample volume from the vial. The syringe plunger moves up to introduce the lower air gap. The lower air gap is within the needle and is not visible.
6. The carousel moves forward so that the syringe needle is over the injector. The syringe is lowered into the injector and the sample expelled.

7. The syringe carriage moves up to withdraw the needle from the injector, and the syringe is washed for the length of time specified in the method. You may select a wash solvent from either reservoir A, reservoir B, or reservoir A and reservoir B (reservoir A then reservoir B). If you choose to wash from reservoir A then reservoir B, the system defaults to reservoir A after the wash cycle is complete. This means that you must fill reservoir A with the solvent you want to use for your solvent plug.
8. If a second sample injection is called for in the sample list, the cycle begins again at Step 1.

NOTE: The wash cycle occurs during the run and does not affect sample turnaround time, except for samples with very short analysis times.

SAMPLE VIALS

The 48-vial carousel takes a standard 2-mL glass screw top vial with outside dimensions of 12 x 32 mm (12 x 35 mm with cap). The vial is sealed with a plastic screw cap with a 5-mm hole, and a white silicone rubber septum (8-mm diameter) faced with an inert film on one side. The inert side is the colored side of the septum. Place the colored side of the septum toward the sample solution.

- **Selecting Sample Vials.** Various sample vials are available from a number of sources. However, the vials you use with the AutoSampler must have the same external dimensions as the vials available from Varian (see Section 13, **8200 Parts and Supplies**).

The volume of sample available is likely the primary consideration when deciding on the appropriate sample vial. When the sample volume is very limited, use vials with very low volume, such as machined glass conical vials (80-100 μ L), or disposable microvial inserts with screw caps or crimp top closures (minimum target diameter of 5 mm) that fit inside a standard 2-mL vial. Refer to the Varian Analytical Supplies Catalog for sample vials that best fit your application.

- **Installing Vial Septa.** The 2-mL vials (P/N 03-949835-00) are shipped complete with plastic screw caps and Teflon®-laminated silicone rubber septa. Always place the colored, chemically inert Teflon laminate face down (next to the contents of the vial).
- **Changing Vial Septa.** The plain glass vial with cap and septum are inexpensive and are often discarded after analysis. You will not need to replace the septa for these vials. The machined glass microvials, however, are relatively expensive and you may want to reuse these vials with new septa (P/N 69-000169-00, 1 gross lot).

Section 5

Performing Injections

INJECTING A SINGLE SAMPLE


Before injecting a sample, be sure you have

- Turned on the AutoSampler. The power switch is located on the rear of the tower.
- Configured the AutoSampler in the Star Chromatography Workstation.
- Built your Method(s) for control of your 8200 AutoSampler and 3800 GC.


ACTIVATE A METHOD

To activate a method, proceed as follows:


1. Turn on the AutoSampler and the GC.
2. From the Start page click on the System Control/Automations button or from the Windows Start button, then click on the System Control icon in the Varian Star 4.5 folder. The System Control Configuration screen appears.




3. Click on the appropriate selection under Instrument in the menu bar for the instrument you wish to inject your sample into.




The System Control Instrument window appears.



4. Click on Method File under File in the menu bar at the top of the System Control Instrument window. Then click on Activate.




5. The Activate a System Control Method File dialog box appears. This dialog box defaults to displaying only files with the .mth extension. Type the name of the method file you wish to activate and click on OK, or scroll through the list of files, and double-click on the name of the method file.




NOTE: You can view the Method Notes by selecting Notes... in the Method File dialog box. If you want the Notes window to appear every time you change and save a Method, check the “Prompt for Notes on Save” box in the Method Notes window.

6. A dialog box asks you if you want to inject or recalculate samples using this method. Click on Inject.




7. Wait for “Ready” to appear in the title bar of the System Control window.




INJECT THE SAMPLE

To inject the sample, proceed as follows:


1. Click on “Inject Single Sample...” under Inject in the menu bar.



2. The Instrument Parameters dialog box appears.




- You may change the Instrument and Operator fields if you wish.
 - The Max Errors field sets the maximum number of non-fatal errors the Instrument will tolerate. The default value is 0.
 - Selecting “Prompt on Sequence Start” will cause this Instrument Parameters dialog box to be displayed whenever you perform an injection or begin a sequence. If you do not wish this window to appear, deselect this item.
3. If the Instrument Parameters dialog box is open, click on **OK**.
 4. The Inject 8200 Sample dialog box appears. Enter the sample name or other identification, along with the number of injections and the vial number you wish to inject from.



If you want to change the Method to be used for the injection, use the Browse button to bring up the **Select a Method File** window.

5. Click on the Inject button. The AutoSampler begins sampling the vial you specified and then injects from it.
 - The Star Chromatography Workstation builds two files once you have clicked Inject: MANUALn.SEQ and MANUALn.SMP. These files replace any previously active Sequence and SampleList files. MANUALn.SMP stores the data you entered in the Inject Sample dialog box.
 - Status messages appear on the message bar at the bottom of the System Control window during the injection.
6. To display the Sequence log when the injection has completed, double-click on the message bar at the bottom of the System Control window.



CREATING AND RUNNING A SAMPLELIST

The SampleList allows you to identify and set different parameters for one to hundreds of samples. The SampleList is intended for use in Automation when performing large numbers of injections. The Star Chromatography Workstation always has an active SampleList associated with it. You activate SampleLists in much the same way you activate Methods. You can also edit a SampleList.


Before building a SampleList, be sure you have

- Turned on the AutoSampler and the GC.
- Configured the AutoSampler to your GC. This is described in **Installing the 8200 Star Workstation Driver** on Page 25.
- Built your Method(s) for control of your 8200 AutoSampler and 3800 GC.


BUILD THE SAMPLELIST

To build a SampleList, proceed as follows:

1. From the Start page click on the System Control/Automations button, or from the Windows Start button click on the System Control icon in the Varian Star 4.5 folder. The System Control Configuration screen appears.
2. Click on the appropriate selection under Instrument in the menu bar for the instrument you wish to inject your sample into. The System Control Instrument window appears.
3. Click on “Active SampleList” under Windows in the menu bar.



4. The Active SampleList appears.
 - The Active SampleList will be named STARTUPn.SMP where n is the number of your instrument quadrant. This SampleList should initially display a single line.



- To resize the SampleList, click on either the borders or corners of the SampleList and drag it to the desired size. You should ordinarily resize your SampleList so that it extends from the left to right sides of your monitor.
- To scroll the cells to the right or to the left, use the scroll bar under the right half of the table.
- To add a line, click on the Add button on the right side of the SampleList. A default entry is appended to the SampleList.

FILL IN THE SAMPLELIST

To fill in the SampleList, proceed as follows:


1. With the Active SampleList window open, type the name or other identification of the sample you wish to inject in the Sample Name field.

	Sample Name	Sample Type	Inj.	Injection Notes	AutoLink	Vial	Injection Volume	Automix Routines
1	Sample-1	Analysis	2	none	none	1	1.0	none
2	Sample-2	Analysis	2	none	none	2	1.0	none
3	Sample-3	Calibration	2	none	none	3	1.0	none
4	Sample-4	Baseline	2	none	none	4	1.0	none
5		Print Calb						
6		New Calb Block						
7								
8								
9								
10								


2. Select sample type: Analysis, Calibration, Verification, or Baseline.
3. Enter the number of injections to be made from the sample in the injection (Inj.) field.
4. Enter information or comments about the sample in the **Injection Notes** field.
5. Verify the vial position on the carrousel that you will be putting your sample into.
6. Enter the volume to be injected in the **Injection Volume** field if you are doing liquid injections.

The Copy Carrousel feature can speed up the process of SampleList building by automatically entering lines for each vial found in the carrousel. The vial position will be filled in for you along with your SampleList defaults.


1. Place a few septum-capped vials in your autosampler carousel and click the Carousel button. The Copy Carousel window appears.



2. Click on the Copy button. The Scanning Carousel window appears.



3. The AutoSampler will rotate and scan the Carousel. When scanning is finished, a line will be entered into the SampleList for each vial position in the Carousel.




HOW TO SET THE DEFAULT SETTINGS


To set the default settings, proceed as follows.

1. With the Active SampleList window open, modify a line of your SampleList to reflect the desired default conditions.
2. Click on the line number associated with the line that you just modified.

3. Click on “Save As Default Sample.”
4. A window asks you whether you wish to save the save SampleList line as the default sample. Click **Yes**.




5. The new lines that you now add will have the same values that you set in the default line. You may change the default settings at any time by repeating this procedure.




HOW TO DELETE SAMPLELIST ENTRIES

1. With the Active SampleList window open, click on the number at the left side of the SampleList beside the first entry you wish to delete.
2. Drag down the number column until you have highlighted all of the entries you wish to delete.
3. Release the mouse button and then click on the Delete button at the right side of the SampleList.




- All of the highlighted lines are removed.




DOWNLOADING A METHOD FROM THE SAMPLE LIST

To download a method and associate it with samples in the SampleList, proceed as follows:


1. Build a SampleList as described in **Build the SampleList** on Page 56.
2. Insert an entry before the desired sample and click on the Sample Type column cell. The cell becomes highlighted.
3. Click on the arrow at the right of the cell to bring up a pull-down menu, then click on Activate Method.




4. Click on the cell under the AutoLink column across from the Sample Type cell selected in Step 2. The Activate Method window opens.



- Select the Method path name to be activated by using the Browse button or typing in the name.



- Click on **OK**. The Method is entered in the SampleList.



- Fill in the SampleList in the lines following for the samples you wish to associate with the method of Step 6.
- Repeat steps 1 through 8 for each series of samples. See the following diagram.


Method Download Lines

	Sample Name	Sample Type	Inj.	Injection Notes	AutoLink	Vial	Injection Volume	Automix Routines
1	Sample - 1	Analysis	5	none	none	1	1.0	none
2	Sample - 2	Analysis	5	none	none	2	1.0	none
3	Sample - 3	Activate Method	5	none	minimum.mth			
4	Sample - 4	Analysis	5	none	none	4	1.0	none
5	Sample - 5	Analysis	5	none	none	5	1.0	none
6	Sample - 6	Activate Method	5	none	111.mth			
7	Sample - 7	Analysis	5	none	none	7	1.0	none
8	Sample - 8	Analysis	5	none	none	8	1.0	none
9	Sample - 9	Activate Method	5	none	startup1.mth			
10	Sample - 10	Analysis	5	none	none	10	1.0	none
11	Sample - 11	Analysis	5	none	none	11	1.0	none

BEGIN THE SAMPLELIST

You are now ready to inject the samples in your SampleList. Since you have been using the Active SampleList for the instrument you have configured, there is no need to save your edits. They are already saved.

1. Click on the Begin button.
2. If the Instrument Parameters dialog box is displayed, click on **OK**. The Begin SampleList window appears.



This dialog box will show the current active Instrument Method. If you wish, you can select a different Method to use at this time by using the Browse... button. When you have selected your Method, click on OK. The 8200 AutoSampler will search for the vials you have specified in your SampleList and make the injections.

BUILDING AND RUNNING A SEQUENCE

A Sequence is essentially a list of Methods and SampleLists that will be run together under automation. The Star Chromatography Workstation always has an active Sequence associated with it. When you activate and subsequently begin a sequence, there is no need for you to individually activate either the Method or the SampleList since the Sequence does that for you.


Before building a Sequence, be sure you have

- Turned on the AutoSampler and the GC.
- Configured the AutoSampler to your GC. This was described on Page 27, **Installing the 8200 Star Workstation Driver**.
- Built your method(s).
- Built the SampleList(s). This was described in **Creating and Running a SampleList** on Page 55.


OPEN THE ACTIVE SEQUENCE WINDOW

To open the Active Sequence window, proceed as follows:

1. From the Start page click on the System Control/Automations button, or from the Windows Start button click on the System Control icon in the Varian Star 4.5 folder. The System Control Configuration screen appears.
2. Click on the appropriate selection under Instrument in the menu bar for the instrument you wish to inject your sample into. The System Control Instrument window appears.
3. Click on “Active Sequence” under Windows in the menu bar at the top of the System Control screen.




- The Active Sequence window appears.



HOW TO ADD A LINE

To add a line to the Sequence, proceed as follows:

1. With the Active Sequence window open, click on the Add button located on the right side of the Sequence.
 - A default entry consisting of the Action Inject, the Method STARTUPn.MTH and the SampleList/Log STARTUPn.SMP appends to the Sequence.



SET THE ACTION, METHOD AND SAMPLELIST/LOG

To set the Action, Method, and SampleList/Log, proceed as follows:

1. Set the Action to "Inject" by clicking on the cell to be modified, scrolling through the list of actions until you come to "Inject," and then clicking on "Inject."
2. Set the Method by clicking once on the Method cell and then clicking on Browse. From the Select a Method File dialog box, if you are not in the Star directory, change directories to Star, then double-click on the name of the method (.mth file) you wish to insert.
3. Set the SampleList/Log in the same way: click on the SampleList/Log cell, and click on Browse. From the Select a Method File dialog box, if you are not in the Star directory, change directories to Star, then double-click on the name of the sample list (.smp file) you wish to insert.

BEGIN YOUR SEQUENCE

To begin your sequence, proceed as follows:

1. When the System Control window indicates Ready at the top of the window, click on the Begin button at the bottom of the active sequence window.

2. Change the fields in the Instrument Parameters dialog box as appropriate. A dialog box tells you the next actions the Star Chromatography Workstation will take.
3. Click on **OK**. A series of messages appears in the message bar at the bottom of the System Control window.
4. Double-click on the message bar. The Sequence log appears, indicating that the selected Method and SampleList have been activated. The Begin button is grayed out and the Suspend button is active.
5. When the Sequence has completed, the header in the System Control window again displays Ready.

AUTOMIX

Automix is a feature available in the SampleList when performing liquid injections. When performing an Automix routine, the AutoSampler extracts a specified amount of liquid from a vial, and then transfers that liquid to the sample vial from which the injection is to be made.

1. The AutoSampler must be configured for a 10- μ L syringe.

NOTE: Do not use Automix with the 100- μ L syringe.

To run Automix, proceed as follows:

1. Build a SampleList as described in **Creating and Running a SampleList** on Page 55.
2. Fill in the SampleList until you reach the Automix Routines column.
3. Click on box under the Automix Routines column for the sample you wish to mix. The box will read 'none.'

	Sample Name	Sample Type	Inj.	Injection Notes	AutoLink	Vial	Injection Volume	Automix Routines
1	Sample - 1	Analysis	5	none	none	1	1.0	none
2	Sample - 2	Analysis	5	none	none	2	1.0	none
3								
4								
5								
6								
7								
8								
9								
10								

4. Click on **none** in the Automix Routines column. The Automix Routines dialog box appears.
5. Click on the first line of the Automix Routines window. The default parameters appear.

	Automix Name	Vial	Volume	Mixes	Mix Stroke	Washes	Reaction Time
1	Add Standards	48	2.0	10	6.0	2	3.5
2							
3							
4							
5							
6							
7							
8							
9							

Buttons: Add, Insert, Delete, Ok, Cancel, Automix First Injection Only, Save As Default Automix

6. Enter a name to describe the function of this Automix routine under Automix Name. You may enter up to 19 characters.
7. Enter the number of the vial number from which the liquid will be extracted under Vial.
8. Enter the amount of liquid to be extracted under Volume. Range: 0.1 to 10.0 μ L.
9. Enter the number of mixes to be performed on the sample vial after liquid has been transferred under Mixes. Each mix withdraws the mixing stroke volume from the sample vial and squirts it back into the vial. Range: 0 to 100.
10. Enter the syringe volume used for mixes under Mix Stroke. Range: 0.1 to 10 μ L.
11. Enter the number of syringe washes to be performed after mixing under Washes. Range: 0 to 9.
12. Enter the delay time following mixing before the sample is injected under Reaction Time. Range: 0.0 to 99.9 minutes.
13. Click on **OK** to return to the SampleList window.

Active 8200 Standalone SampleList: STARTUP1.SMP

	Sample Name	Sample Type	Inj.	Injection Notes	AutoLink	Vial	Injection Volume	Automix Routines	
1	Sample - 1	Analysis	5	none	none	1	1.0	1 routine	Add
2	Sample - 2	Analysis	5	none	none	2	1.0	none	Insert
3									Delete
4									
5									
6									
7									
8									
9									
10									

Begin Suspend Resume Carrousel... Save As Default Record

When this sample is processed, the following takes place:

1. The AutoSampler selects the vial specified in the Automix Routines window.
2. The syringe withdraws a volume equal to the amount specified in the Automix Routines window from the Automix Routines vial and transfers the contents to the vial specified in the SampleList window.
3. The syringe mixes the sample by withdrawing the amount specified in the Automix Routines window from the SampleList vial and squirting it back into that vial. This is done the number of times specified in the Automix Routines window.
4. The syringe is washed the number of times specified in the Automix Routines window. Countdown of the reaction time specified in the Automix Routines window begins.
5. After the reaction time countdown has completed from the last Automix routine, the syringe withdraws and injects the volume specified in the SampleList window from the SampleList vial into the GC.

Section 6

SPME Injection Setup

To set up the 8200 AutoSampler prior to using the SPME Injection Method, you will

- Build a SPME Injection Method
- Install the fiber holder on the syringe carriage.
- Align the fiber with the injector port.
- Adjust the vial retainer arm if you will be using the 12-vial carousel.

None of these individual tasks should take you more than five minutes to complete.

BUILD A NEW SPME INJECTION METHOD


NOTE: Be sure you have selected the **SPME Device Setting** in the **8200 Hardware Configuration window (Section 2, Specify the 8200 Hardware Configuration)**.

The method is a record of all the instructions and parameters required to analyze your sample. The Star Chromatography Workstation allows you to create separate methods for each of your sample types.

To build a new SPME Injection Method, proceed as follows:


1. From the Start page click on the Method Editor button, or from the Windows Start button click on the Method Editor icon in the Varian Star 4.5 folder.

The Method Editor window appears with the Create/Open Method File in the foreground.



2. Select "Create a new method file" and click **OK**.

- The New Method Section window appears. Select “8200 Stand-Alone AutoSampler” in the Instrument Control section of the window.



- Referring to the following table, use the up and down arrows to set the Bus Address. Click on **OK**.

If the 3800 GC Address is...	Set the Bus Address
44	26
45	30
46	34
47	38

- The New Method Section window closes and the 8200 Configuration icon appears in the lower left hand corner of the Method Editor window.
- Double-click on the 8200 Configuration icon. The Configuration window appears.

SPECIFY THE CONFIGURATION WINDOW OPTIONS

- Select the SPME Mode option under Sampling Modes in the Configuration window.

Configuration

Carrousel Type: 48 Vials

AutoSampler Timing

GC Cycle Time: 1.00

Prep Ahead

Sampling Modes

Liquid

Standard Mode Viscous Sample

Volatile Sample User Defined

Neat Sample

SPME Mode

SPME Options

Absorb Time (min): 10.00

Desorb Time (min): 2.00

Sample Headspace Agitate


Save Cancel

NOTE: The SPME III Agitation option is an accessory option which provides for agitation of the sample during absorption time, thereby increasing the rate at which solutes are extracted by the SPME fiber coating. See Section 8, Operation of the SPME Agitation Accessory.

2. Set the Carrousel Type to 48 Vials or 12 Vials using the pull-down menu that appears when you click on the arrow to the right of the box. The default value is 48 Vials.
3. Determine the length of a GC cycle (run time plus cooldown plus re-equilibration) by allowing the GC to cycle two or three times.
4. Under AutoSampler Timing, set the GC Cycle Time by entering the time in minutes.
5. Select Prep Ahead. Use of this feature is highly recommended with SPME to minimize dead time between runs. It allows absorption onto the SPME fiber of a sample while the previous sample is undergoing a chromatographic run.
6. Under SPME Options, set the Absorb Time in minutes.
7. Under SPME Options, set the Desorb Time in minutes.
8. If you are going to use the SPME Injection Method to sample headspace, check Sample Headspace by clicking on the box.
9. Click on **Save**. The Configuration window closes.

SAVE YOUR METHOD

1. Close the Method Editor window. A window appears that asks you whether you want to save the changes you have made to the method.
2. Click **Yes**. The Save Method File As window appears prompting you to name the new method.



3. Assign a name of up to eight letters and/or numbers to the method, and click **OK**. The .mth extension will be added automatically.

INSTALL THE FIBER HOLDER ON THE SYRINGE CARRIAGE

To install the fiber holder on the syringe carriage, proceed as follows:

1. Remove the backing from the adhesive part of the AutoSampler's solvent line clamp (P/N 22-120202-00) and press the clamp into place (one-inch directly above the solvent line hole on the carriage).
2. Open the front door of the syringe mounting assembly by *simultaneously* and firmly lifting *both* side locking tabs. Raise the tabs to the upper position. This may require considerable force, so hold the carriage assembly down with one hand. Alternatively, you may hold the mount down with both thumbs while you lift the tabs. The locking mechanism will become less stiff with use. A red tab is visible above the syringe mount assembly when the door is unlocked.

- Remove the liquid syringe and place the solvent line in the solvent line clamp.




**WARNING:
EYE HAZARD**


Wear Proper Eye Protection. Solvent may spray out of the Teflon® tubing when disconnected from the syringe.

- Place the fiber into the fiber holder according to the installation procedures supplied with Supelco™ SPME fiber and fiber holder. Refer to the following diagrams.

SPME Holder with Fiber Assembly



Fiber Assembly



- A SPME Fiber
- B Septum Piercing Needle
- C Needle Ferrule
- D Retaining Nut
- E Fiber Attachment Needle
- F Barrel
- G Plunger
- H Sealing Septum
- I Color-Coded Screw Hub
- J Slot
- K Retaining Screw

NOTE: If a fiber is already installed in the fiber holder, make sure that the plunger is withdrawn so that the fiber is completely inside the septum piercing needle before disassembling any part of the SPME holder.

- Install the fiber assembly in the AutoSampler syringe mount.

NOTE: If the agitation accessory is installed, see Section 8 for details on fiber holder installation.

NOTE: If, after you have installed the fiber assembly, the fiber extends from the sheath, do not proceed with the fiber alignment until you have first activated your SPME injection method. To do this, follow the steps given on Page 51, Activate a Method, using the method you created in Build a New SPME Injection Method, Page 69. This will cause the fiber to retract into the sheath even if you have not yet completed the remaining steps in this SPME injection method setup.

6. Close the syringe mount door.
 - Check that the door is fully closed before pushing the side-locking tabs down. Then, holding the door closed and supporting the carriage with one hand, press the side locking tabs down with the other hand. This locks the syringe into position so that the syringe hex screw does not move. When the tabs are down, the door is locked and the red signal tab is not visible.

ALIGN THE FIBER WITH THE INJECTOR PORT

When you initially install the fiber in the syringe mount, the alignment of the fiber with the needle port of the injector nut may not be precise.


NOTE: Once the fiber is correctly mounted, you will not need to adjust it again.

To align the fiber with the injector port, proceed as follows:

NOTE: Be sure that the injector insert is compatible with the SPME injection technique, i.e., do not use a packed insert or high performance insert.

NOTE: The GC should be cooled down during this operation.

1. Check that the AutoSampler power switch is OFF. The power switch is located on the rear of the tower.
2. Press down on the storage module latch release and swing the module out to access the GC injector. When the storage module swung out (with power off), a contact in the latch assembly activates that turns the syringe carriage stop solenoid on or off each time the carriage position sensor is interrupted.



NOTE: You must press the latch release down when relatching the storage module.

3. Wearing gloves to prevent contamination, remove the injector nut and replace the glass insert in the injector, as described in your injector manual.
4. Replace the injector nut, and verify that the column head pressure is adequate.
5. Turn the GC power ON.
6. Turn the AutoSampler power ON.
7. Raise the syringe carriage from its rest position to its top limit. The stop solenoid toggles (you will hear a “click”) to allow you to lower the carriage by hand.



Take care that the carriage does not drop and damage the fiber assembly.

8. Loosen the screws holding the mounting plates in place so that the AutoSampler can be moved on the GC cover.
9. Lower the syringe carriage by hand until the fiber/sheath tip is just above the GC injector needle port.
10. Carefully adjust the position of the AutoSampler to bring the fiber/sheath tip into precise alignment with the injector needle port.
11. Securely tighten the screws holding the mounting plates in place.

12. Lower the syringe carriage and puncture the injector septum three or four times with the fiber sheath to ensure that the fiber is properly aligned with the injector port.
 - Repeat this procedure each time you replace the injector septum.
13. To return the syringe carriage to its rest position, raise the syringe carriage to its upper limit of travel. Then lower the carriage a fraction of an inch until it engages.
14. Close the storage module.

NOTE: To re-latch the storage module, the latch release must be pressed down.

SPME FIBER CALIBRATION PROCEDURE

The Varian 8200 SPME AutoSampler was designed to minimize the need to make adjustments to the hardware after the initial installation. However, some users have noticed that the position of the SPME fiber within the protective fiber sheath is not optimum after a fiber has been replaced.

To achieve the best performance after a new fiber has been installed on an SPME autosampler, the position of the fiber within the sheath should be adjusted. The fiber in its rest position should be neither protruding from the sheath, nor be too far up inside the sheath. Ideally, the tip of the fiber will be even with the end of the sheath. To perform this adjustment, proceed as follows:

First, determine the module number used by the 8200. When the 8200 Instrument window is displayed on your screen, the module number is contained in the header as the suffix; for example, the header "8200.28" identifies this 8200 as module 28. The module number can also be found in the icon label. All changes described below must be done for the appropriate module.

1. Make certain that the SPME syringe is specified as Hardware in the 8200 System Control page. Exit System Control. (Do not simply minimize the program, but exit it altogether. System parameters are only read at the time of program startup.)
2. Find the file named WS8200.INI in the WINDOWS directory. Use "Notepad" to edit this file.
3. Find the appropriate section for your instrument in the WS8200.INI file. For example, if your 8200 AutoSampler is module 28, all parameters affecting this module will be found under the header [8200.28]. Add the following line to this section:

PlungerSpmeRetracted=150

4. Save and exit WS8200.INI
5. Startup System Control. As the autosampler initializes itself, you will see that the fiber is now extended from the sheath. Measure this extended portion of fiber with a millimeter scale.
6. Multiply the length of extended fiber by 7.56. For example, if the fiber measurement were 4.8 mm, then $4.8 \times 7.56 = 36.3$. Round off the result to the nearest whole number (36 in the example), and add it to 150. The example would give a value of 186.
7. Exit System Control. Edit the PlungerSpmeRetracted parameter in the WS8200.INI file to change the value(s) to the number just calculated. In the example, this value was 186. Save the file(s) and exit.
8. Startup System Control. After the autosampler initializes, the fiber should be in the correct position within the sheath.
9. If the position of the fiber is still not correct, change the value of the PlungerSpmeRetracted parameter incrementally by 4 units at a time until the position is optimum. A change of 4 units will move the fiber up or down approximately 0.5 mm.

ADJUST THE VIAL RETAINER ARM HEIGHT


To adjust the vial retainer arm height, proceed as follows.

1. Locate the vial retainer arm on the AutoSampler.

- A Solvent Line
- B Plunger Drive
- C Carriage
- D Side Locking Tabs (2 places)
(Down: door closed; red tab not visible)
- E Syringe Mount
- F Carriage Bracket
- G Needle Guide
- H Vial Retainer Arm




NOTE: Be sure this label is placed in full view on lower portion of carriage.



2. Using a 0.05-inch hex wrench, loosen the hex screw that holds the vial retainer arm in position.
3. Slide the vial retainer arm up or down to the stop clip to accommodate 10-mL or 2-mL vials, respectively. Check that vial clears retainer arm and readjust for slight clearance.
4. Tighten the hex screw.
5. Check arm alignment on Page 25.

- A Vial Retainer Arm Block
- B Hex Screw
- C 10 mL Vials (adjust up)
2 mL Vials (adjust down)
- D Vial Retainer Arm



FINAL SPME INJECTION METHOD SETUP CHECKLIST

Setup of the SPME Injection Method for your 8200 AutoSampler is now complete. Before operating the AutoSampler as described in Section 7, make sure you have

- Built a SPME Instrument Control Method
- Installed the fiber on the syringe carriage
- Aligned the fiber with the injector port
- Adjusted the vial retainer arm for the size of vial you will use
- Installed a SPME compatible injector insert

SPME ADJUSTMENTS

The Varian 8200 SPME AutoSampler was designed to minimize the need to make adjustments to the hardware after the initial installation. However, some users need to “fine-tune” the system, either because of their special sampling requirements, or because of particular hardware requirements. This document describes a procedure which can be used to change a number of the sampling parameters.



Changing the parameters described in this document can lead to damage to your autosampler. The 8200 AutoSampler is a robotic system which is programmed by the following parameters. The standard software supplied with the instrument has been thoroughly tested to ensure that no conflicts occur during operation. However, if the preset parameters are altered, then problems can occur. For example, changing the distance the fiber is inserted into the vial can break the fiber if it hits the bottom of the vial. For this reason, make small changes only. Larger changes should only be made incrementally, confirming correct operation between each change.

One file controls all the sampling parameters: WS8200.INI. This file is found in the WINDOWS directory, usually located on your C: hard drive. It can be edited with any text editor or word processor. However, if a word processor is used, be sure to save the file as plain text, NOT as a *.DOC file. In addition, the name of the file must not be changed in any way.

The parameters that control SPME sampling are:

- PlungerSpmeSensorCount=22
- PlungerSpmeRetracted=200
- PlungerSpmeStepsIntoVial=100
- PlungerSpmeStepsIntoInjector=100
- PlungerSpmeSpeed=120
- SpmeHeadSpaceDepthPct=5
- SpmeLiquidDepthPct=36
- SpmeIntoInjector=730

To make changes to parameters, use the following procedure:

First, determine the module number used by the 8200. When the 8200 Instrument window is displayed on your screen, the module number is contained in the header as the suffix; for example, the header "8200.28" identifies this 8200 as module 28. The module number can also be found in the icon label. All changes described below must be done for the appropriate module.

1. Make certain that the SPME syringe is specified as Hardware in the 8200 System Control page. Exit System Control. (Do not simply minimize the program, but exit it altogether. System parameters are only read at the time of program startup.)
2. Find the file named WS8200.INI in the WINDOWS directory. Use "Notepad" to edit this file.
3. Find the appropriate section for your instrument in the WS8200.INI file. For example, if your 8200 AutoSampler is module 28, all parameters affecting this module will be found under the header [8200.28]. Add or modify lines in this section. For example, you might want to increase the sampling depth of the fiber in the vial. The parameter that controls this is:

PlungerSpmeStepsIntoVial=100

The default is 100; increasing this number will cause the fiber to be pushed further down into the sample after the fiber sheath has penetrated the vial septum. Remember to change this value by small amounts; for example, changing the value to 108 will change the insertion depth by about 1 millimeter. See the PARAMETERS section below for more information on each of the sampling parameters.

4. Save and exit WS8200.INI
5. Startup System Control.
6. Start up the AutoSampler and observe whether the changes you made were adequate for your purpose. If all is correct, stop here; otherwise, continue to the next step.
7. Exit System Control. Re-edit the appropriate parameters in the WS8200.INI file to change the value(s) as needed. Save the file(s), exit, and go back to Step 5.

SPME PARAMETERS

PlungerSpmeSensorCount. This parameter is used to set a zero position for the SPME syringe plunger. It should normally not be modified.

PlungerSpmeRetracted. Used to determine the rest position of the SPME fiber within the protective sheath. This parameter is the one most frequently modified by the user, because small differences between autosamplers, SPME holders, and fibers can affect the rest position of the fiber. To set this parameter correctly, see the procedure for calibrating the fiber position on Page 76.

PlungerSpmeStepsIntoVial. This parameter determines how far down the SPME fiber is extended from the sheath after the sheath has penetrated the vial septum. Increasing this number will extend the fiber further. A too-large value will result in the fiber breaking on the bottom of the vial.

PlungerSpmeStepsIntoInjector. Used to determine how far the SPME fiber will be extended after the carriage has moved down to the injector and the fiber sheath has penetrated the injector septum. The higher the number, the further the fiber will be extended. This parameter can be used to adjust the desorption position of the fiber within the injector.

PlungerSpmeSpeed. The rate at which the fiber plunger moves up or down. This number should normally never need to be changed. A value of 120 corresponds to a rate of approximately 16 mm per second.

SpmeHeadspaceDepthPct. This parameter determines how far into the vial the sheath will be inserted before the SPME fiber is extended for adsorption during headspace sampling. The default value of 5 allows the fiber sheath to pass just through the vial septum before fiber extension occurs. Use of certain vials with slightly different dimensions may lead to fiber breakage, because the sheath does not penetrate the vial septum. This can be corrected by increasing the value of this parameter slightly. The range is 0 to 100, with 0 representing the top of the vial.

SpmeLiquidDepthPct. Determines how far into the vial the sheath is inserted before fiber extension during SPME liquid sampling. The range of values is exactly the same as for the parameter, SpmeHeadspaceDepthPct. Using a value of 100 will cause the sheath to go all the way to the bottom of the vial. In this case, the fiber would be broken upon extension.

SpmeIntoInjector. Determines how far the SPME sheath is inserted into the injector before the fiber is extended for desorption. This would be modified only if injectors have a non-standard height. The value of this parameter was determined for injectors whose septum nuts are even with the GC cover. If an injector is below this point, for example, then the value of this parameter should be increased.

Section 7

SPME Injection Operation

This section gives additional information that you will need to operate the 8200 AutoSampler using the SPME Injection Method. After familiarizing yourself with the information in this section, proceed directly to Section 5, **Performing Injections**.

SAMPLING

LOADING THE CARROUSEL

The carousel rotates clockwise and the vial position numbers increase counter-clockwise around the carousel. Observe this order when loading sample vials. Load sample vials according to the vial numbers shown on the carousel label.


The AutoSampler ignores vials without a septum. Seal all vials in the carousel that are to be sampled with a cap and septum.

SAMPLE VIALS

The 48-vial carousel takes a standard 2-mL glass Screw Top vial with outside dimensions of 12 x 32 mm (12 x 35 mm with cap). The vial is sealed with a plastic screw cap with a 5-mm hole, and a white silicone rubber septum (8-mm diameter) faced with an inert film on one side. The inert side is the colored side of the septum. Place the colored side of the septum toward the sample solution.

NOTE: The 48-vial tray can be placed on the turntable in any position.

The 12-vial carousel takes a standard 10-mL glass Crimp Top vial with outside dimensions of 25 x 50 mm (25 x 54 mm with cap). The vial is sealed with a crimp top cap with a 9.5-mm hole, a Viton® septum, and a silicone-rubber o-ring. It is recommended that you place the o-ring in contact with the crimp top, and the Viton® septa between the o-ring and the top of the vial.



NOTE: Other vials compatible with the Varian 8200 12-vial tray will become available later. Check with your Varian representative for availability.

NOTE: The 12-vial tray must be placed on the turntable with the pin aligned with vial 1.

- **Selecting Sample Vials.** The vials you use with the AutoSampler must have the same external dimensions as the vials available from Varian (see AutoSampler Parts and Supplies).
- **Installing Vial Septa.** The 2-mL vials (P/N 03-949835-00) are shipped complete with plastic screw caps and Teflon®-laminated silicone rubber septa. Always place the colored, chemically inert Teflon laminate face down (next to the contents of the vial).

The 10-mL vials (P/N 03-918873-02) are shipped with crimp top caps with black Viton® septa inserted. Remove the black Viton septum; place the o-ring in the cap, then replace the thin Viton septum. Place the cap on the vial.

NOTE: To ensure that the SPME fiber experiences a long lifetime, use only septa supplied by Varian.

CONDITIONING THE FIBER

A new fiber must be conditioned by desorbing it for a minimum of 15-20 minutes in an injector that is at least 10°C warmer than the temperature to be used during the analysis. The GC column should then be temperature programmed. This procedure should be repeated until there are no extraneous peaks.

SAMPLE VOLUMES

- **Liquid Sampling.** Fill the 2-mL vials with no more than 1.2 mL of liquid.
Fill the 10-mL vials with exactly 10 mL of liquid. These vials actually hold 16 mL, so a large air space remains over a 10-mL sample.
- **Headspace Sampling.** Fill the 2-mL vials with no more than 0.8 mL of sample.
Fill the 10-mL vials (they actually hold 16 mL) with no more than 10 mL of sample.

SETTING THE SPME/GC PARAMETERS

- **Injector Temperature.** The injector temperature should be the lower of the following: 20°C below the column limit and the SPME fiber temperature maximum. One to two minutes is usually sufficient for desorption, but this value varies with fiber and sample type.
- The recommended injector inserts for SPME are shown below:

Injector	Insert	Part Number
1079	1079 SPME Insert	03-925330-00

RUNNING THE TEST SAMPLE

After installing the SPME software, hardware, and fiber, test the SPME system with the **SPME Sensitivity Test Sample**.

The test sample is composed of 1 ng/ μ L each of nitrobenzene and nitrotoluene in water (1% methanol has been added to stabilize the sample). These compounds were selected because they exhibit a good response with many GC detectors, including the flame ionization detector, the electron capture detector, the thermionic selective detector, and the ion trap detector.

GC Conditions

To run the SPME Sensitivity Test Sample, use the following GC conditions:

Column:	Nearly any capillary column can be used to separate the SPME Sensitivity Test Sample.
Column Conditions:	50°C for 1 minute; then 20°C/minute to 150°C; hold for 2 minutes.
Injector:	200°C isothermal.
Detector:	Settings depend on the detector used.

SPME Software Settings

Adsorb Time: 10 minutes
Desorb Time: 1 minute

SPME CONDITIONS

Fiber: 100 μ m PDMS
Glass Insert: SPME
Carrier Gas: 4.7 mL/min

The test sample is composed of 1 ng/ μ L each of nitrobenzene and nitrotoluene in water (1% methanol has been added to stabilize the sample). These compounds were selected because they exhibit a good response with many GC detectors, including the flame ionization detector, the electron capture detector, the thermionic selective detector, and the ion trap detector.

A representative SPME test sample chromatogram is shown.

GC CONDITIONS

Injector: 200°C


Column: 50°C/min, hold 1 min

Program: 20°C to 110°C, hold 2 min

Detector: 240°C

The chromatogram shown was obtained with an FID:
range 10, attenuation 128.

- 1 Nitrobenzene
- 2 Nitrotoluene



SPME Test Sample Chromatogram

Section 8

SPME III with Agitation


INTRODUCTION

Solid Phase MicroExtraction (SPME) may be used for extraction and concentration of both volatiles and semi-volatile components. Typically, extraction of volatiles is performed with the SPME fiber inserted into the headspace of the sample. Because of the high diffusion rates of volatiles in the gas phase, no agitation of the sample is usually required to achieve good results in a reasonable amount of time.

Semi-volatile components must be extracted with the fiber inserted into the bulk sample. In this case, sampling times can be quite long before equilibrium occurs of the sample with the fiber coating. The higher the molecular weight of the sample, the lower the diffusion rate, and the longer time that is required for sample uptake by the fiber.

If the sample is well-mixed during the SPME process, then diffusion rates in the liquid sample are no longer a factor in sampling time. For this reason, an agitation accessory was developed as an option for the Varian 8200 AutoSampler. The increase in uptake rate for a specific compound, fluoranthene, can be seen below.

Comparison of Fluoranthene Uptake With and Without Agitation



OPERATION OF THE SPME AGITATION ACCESSORY

PRINCIPLES OF OPERATION

The SPME agitation accessory for the Varian 8200 AutoSampler is used during automated Solid Phase MicroExtraction (SPME) of liquid samples. This device provides for agitation of the sample during extraction, thereby increasing the rate at which solutes are absorbed by the SPME fiber coating. Agitation of the sample during SPME allows the solutes to reach equilibrium with the fiber much more rapidly than when the sample is not agitated.

The SPME agitation accessory works by vibrating the SPME fiber protective sheath. Vibration of the sheath results in vibration of the fiber in the liquid sample, thereby causing the fiber to contact fresh (undepleted) sample continuously.


This mechanism of agitation gives better results than either ultrasonic vibration or magnetic stirring using a Teflon®-coated stirbar. Ultrasonic vibration gives efficient mixing, but heats the sample. Magnetic stirbars do not provide the agitation required to mix the liquid layer near the fiber with the rest of the bulk solution.

DESCRIPTION OF THE AGITATION DEVICE

The SPME agitation accessory consists of a small precision motor which drives a small eccentric wheel (cam). The cam is coupled to a rigid arm, which vibrates back and forth when the cam is spinning. This vibration is transferred to the SPME fiber sheath through a silicone rubber septum, which is mounted in the rigid arm. The fiber sheath passes through the center of this septum.

The septum provides a degree of damping to the vibration, so that the sheath is never forced to move far enough to damage either the sheath or the fiber. The amplitude of vibration and the speed of rotation of the motor have been carefully chosen to provide optimum results, and cannot be changed.

The following picture shows the SPME agitation accessory. The agitation system is mounted on the AutoSampler carriage to the right of the syringe mount. Power is supplied to the motor by a cable attached to the motor with a slip-on connector. The entire mechanism is held onto the mounting plate with three screws, which can be loosened for easy removal of the entire agitation mechanism.



- A Power Connector
- B Agitation Motor
- C Cam
- D Septum Holder



INSTALLATION

The SPME III agitation accessory is factory installed. Field upgrades of the agitation accessory may be performed on 8200 AutoSamplers purchased after June 1996.

Installation of the AutoSampler should proceed exactly as described in Section 2 of this manual.

If the 48-vial carousel is to be used, no further modifications are necessary. If the 12-vial carousel is to be used, the AutoSampler must be modified as follows:

1. Change the position of the vial retainer arm to its uppermost setting (see Section 6, Page 77).
2. Attach the vial retainer arm activator strip to the front, left hand side of the AutoSampler carriage. This strip clips under the carriage at its lower end and is held onto the carriage with a single screw.



A Activator Strip

SETUP AND OPERATION

The following procedure must be followed to set up the SPME agitation accessory for sampling and injections:

1. Install a pre-drilled septum into the rigid agitator arm, if one has not already been installed.
2. Install the SPME fiber holder with fiber assembly into the 8200 syringe holder, making certain that the fiber sheath is inserted through the center of the septum in the agitator arm and the needle guide of the syringe holder.
3. Build a method for SPME sampling with the **Agitate** box checked.

SPME sampling will then occur with vibration of the fiber from the time the fiber is extended from the sheath into the sample to the time just before the fiber is withdrawn back into the sheath. Details of the above procedure are given:


INSTALLATION OF AGITATOR SEPTUM

To install a new agitator septum, first remove the SPME fiber holder, if one is installed. Loosen the Phillips head screw holding the septum retainer to the rigid arm. Rotate the retainer enough so that the old septum may be removed, and insert a new pre-drilled septum. Rotate the septum retainer back to its original position and tighten the retaining screw.

NOTE: If a pre-drilled septum is not available, prepare a new septum by pushing an old fiber sheath or similarly sized needle through the center of the septum to create a guide hole. It is very important to locate the center of the septum in this manner before installing the new septum into the rigid arm.

INSTALLATION OF SPME FIBER HOLDER

Open the syringe mounting assembly by lifting the side locking tabs (see Section 6 for details of this operation). Because of the limited access to the right-hand tab, it may be necessary to use a tool, such as a flat-blade screwdriver, to facilitate this operation. Remove the liquid sampling syringe, if one is installed, and clip the solvent line into the solvent line clamp. After first making certain that the fiber is completely withdrawn into the fiber sheath, insert the sheath through the hole in the center of the pre-drilled septum.



A Pre-drilled Septum Hole

Carefully guide the fiber sheath through the lower needle guide while moving the entire holder assembly into the syringe mounting assembly. Ensure that the end of the plunger is inserted into its guide. Close and lock the syringe mount door.

BUILDING A METHOD THAT INCLUDES AGITATION

Refer to Section 6 for details on building a SPME method. To activate agitation, open the Method Editor for the 8200 Stand-Alone AutoSampler Module. In the SPME Options area of the Configuration window a check box labeled **Agitate** will be found. This box must be checked to activate agitation of samples during SPME operation.

The screenshot shows the 'Configuration' dialog box with the following settings:


- Carrousel Type: 48 Vials
- AutoSampler Timing
 - GC Cycle Time: 1.00
 - Prep Ahead:
- Sampling Modes
 - Liquid
 - Standard Mode:
 - Viscous Sample:
 - Volatile Sample:
 - User Defined:
 - Neat Sample:
 - SPME Mode:
- SPME Options
 - Absorb Time (min): 10.00
 - Desorb Time (min): 2.00
 - Sample Headspace:
 - Agitate:

Buttons: Save, Cancel


CHANGING FROM SPME III TO LIQUID SAMPLING

Liquid sampling with the standard 10 μ L syringe (or static headspace or LVI sampling) requires removal of the agitation accessory. Use the following procedure to remove the accessory:

1. Open the syringe holder and remove the SPME fiber holder. Close and lock the syringe holder.
2. Unplug the power cable from the top of the motor.



3. Using a small Phillips head screwdriver, insert screwdriver in the hole in the side of the tower casting. Loosen this screw a couple of turns.



4. Unlatch the carousel module and swing the module away from the AutoSampler tower. Turn the AutoSampler power off, then back on. Lift the AutoSampler carriage upward until an audible click can be heard. The latch that prevents the carriage from dropping has now been pulled back, and the carriage can be lowered to any position.

5. Lower the carriage until the top screw is aligned with the hole in the side of the tower casting. Loosen this screw a couple of turns.



6. Raise the carriage and loosen the bottom screw.
7. The agitation assembly may now be removed from the mounting bracket by pulling straight out.



8. If the 12-vial carousel has been in use, remove the vial retainer arm activator strip and reset the vial retainer arm to its lowest position.

Converting to SPME sampling from liquid sampling is the reverse of the above procedure.

TROUBLESHOOTING AND MAINTENANCE

The SPME agitation accessory has been designed so that it needs little or no maintenance during routine operation, other than changing the agitator damping septum regularly. It is recommended that the septum be replaced at least once every 100 injections. The septum should also be replaced if it begins to look worn or fragmented.

If the agitator does not function when it should, first check your method for proper setup. Remember, the **Agitate** box must be checked and the Method must be active for agitation to occur. If the agitator still does not work, check all cables and electrical connections. If the motor still does not run, then a qualified Varian Customer Support Representative must investigate the problem.

METHOD DEVELOPMENT HINTS

Agitation of the sample during Solid Phase MicroExtraction increases the rate of absorption of analytes by the fiber. Agitation does not increase the total amount of analyte absorbed at equilibrium. In many cases, it is not necessary to reach equilibrium to achieve good analytical results. If the sampling time is consistent from one sample to the next, reasonable peak area reproducibility can be obtained. This should be considered for analytes that take a long time to reach equilibrium with the fiber, even with agitation.

In general, the higher the molecular weight of a solute, the longer it will take to reach equilibrium with the fiber, due to slower diffusion rates into the fiber. The type of fiber chosen will also affect time to equilibrium.

Methods development for a SPME method incorporating agitation primarily requires determination of the sampling time which is sufficient to provide adequate detection limits for solutes with reasonable reproducibility. Since SPME is a sample preparation technique, and not a direct injection of a sample solution, reproducibilities should be comparable or better than other sample preparation techniques, such as solid phase extraction and Soxhlet extraction.

Reproducibilities for such methods are typically in the range of 5% to 25% RSD (relative standard deviation).

Automated methods development can be performed with the 8200 to determine optimum absorption and desorption times by setting up several methods containing the different parameters. These methods can then be run in sequence, and the results analyzed to determine the best parameters.

Because SPME is a competitive equilibrium technique, usually with aqueous samples, poor reproducibilities may result from interactions not encountered in samples dissolved in organic solvents. For example, the glass walls of vials may adsorb components from the sample, introducing error. It may be necessary to add a displacing agent, such as acid, to the sample to eliminate this effect. Keep in mind that SPME samples through chemical interactions, rather than pure physical measurements, when developing methods for your analytes.

Section 9

Ambient Headspace Injection Setup

To set up the 8200 AutoSampler prior to using Ambient Headspace, you will

- Build the Ambient Headspace Injection Method
- Remove any solvent from the solvent lines
- Connect the pneumatics
- Install the 100- μ L headspace syringe
- Align the syringe needle with the injector port

With the exception of removing the solvent from the solvent lines, which may take 15 minutes or longer, none of these tasks should take you more than 5 minutes to complete.

NOTE: Ambient Headspace Sampling requires removal of the SPME agitation accessory. If this device is installed on your AutoSampler, refer to Section 8, Changing From SPME III to Liquid Sampling on Page 94.

BUILD THE AMBIENT HEADSPACE INJECTION METHOD


The method is a record of all the instructions and parameters required to analyze your sample. The Star Chromatography Workstation allows you to create separate methods for each of your sample types.

NOTE: Be sure you have selected the 100 μ L Syringe Setting. See Specify the 8200 Hardware Configuration on Page 27.


To build a new ambient headspace injection method, proceed as follows:

1. From the Start page click on the Method Editor button, or from the Windows Start button click on the Method Editor icon in the Varian 4.5 folder.

The Method Editor window appears with the Create/Open Method File in the foreground.



2. Select “Create a new method file” and click OK.
3. The New Method Section window appears. Select “8200 Stand-Alone AutoSampler” in the Instrument Control section of the window.



4. Referring to the following table, use the up and down arrows to set the Bus Address. Click **OK**.

If the 3800 GC Address is...	Set the Bus Address
44	26
45	30
46	34
47	38

5. The New Method Section window closes and the 8200 Configuration icon appears in the lower left hand corner of the Method Editor window.
6. Double-click on the 8200 Configuration icon. The Configuration window appears.

SPECIFY THE CONFIGURATION WINDOW OPTIONS

1. Select the User Defined option under Liquid in the Sampling Modes section of the Configuration window.

The screenshot shows the 'Configuration' dialog box with the following settings:

- Carrousel Type:** 48 Vials
- AutoSampler Timing:** End Time: 1.00; Prep Ahead:
- Sampling Modes:** Liquid: Standard Mode, Viscous Sample, Volatile Sample, User Defined, Neat Sample, SPME Mode
- Liquid Sampling Options:** Solvent A Wash; 100 uL Syringe
- User Defined Parameters:**
 - Solvent Flush Sampling: Yes, No
 - Syringe Wash Time (sec): 20
 - Air Dry After Wash: Yes, No
 - Solvent Plug Size (ul): 10
 - Vial Needle Depth (%): 90
 - Uptake Speed (ul/sec): 50
 - Upper Air Gap: Yes, No
 - Lower Air Gap: Yes, No
 - Pause Time (sec): 2
 - Hot Needle Time (min): 0.00
 - Injection Rate (ul/sec): 50
 - Needle Residence Time (min): 0.00


2. Set the Carrousel Type to 48 Vials using the pull-down menu that appears when you click on the arrow to the right of the box.
3. Under AutoSampler Timing, set the End Time using the up and down arrows.
4. The Prep Ahead option should be deselected for Ambient Headspace injections.
5. Under Liquid Sampling Options, set the top selection to “Solvent A Wash” if the solvent A reservoir contains only air, to “Solvent B Wash” if the solvent B reservoir contains only air, or to “A then B” if both reservoirs contain only air. You select this option using the pull-down menu to the right of the box.
6. Under Liquid Sampling Options, set the bottom selection to “100 μ L Syringe” using the pull-down menu that appears when you click on the arrow to the right of the box.
7. Under User Defined Parameters, select “Yes” for Solvent Flush Sampling.
8. Under User Defined Parameters, set Syringe Wash Time (sec) using the up and down arrows.
9. Under User Defined Parameters, select “No” for Air Dry After Wash.
10. Under User Defined Parameters, set the Vial Needle Depth by clicking on the arrows to change the default value of 90%. You should set the vial needle depth so that the needle does not penetrate the liquid in the sample vile. For a volume of 1 mL, a vial needle depth of 25% is reasonable.

11. Under User Defined Parameters, set the Uptake Speed ($\mu\text{L}/\text{sec}$) using the up and down arrows.
12. Under User Defined Parameters, set the Pause Time (sec) using the up and down arrows.
13. Under User Defined Parameters, set the Hot Needle Time (sec) using the up and down arrows.
14. Under User Defined Parameters, set the Injection Rate ($\mu\text{L}/\text{sec}$) using the up and down arrows.
15. Under User Defined Parameters, set the Needle Residence Time (min) using the up and down arrows.
16. Click on **Save**. The Configuration window closes.

SAVE YOUR METHOD

1. Close the Method Editor window. A window appears that asks you whether you want to save the changes you have made to the method.

Click **Yes**. The Save Method File As window appears prompting you to name the new method.



2. Assign a name of up to eight letters and/or numbers to the method, and click **OK**. The .mth extension will be added automatically.


If...	Then....
this is a new AutoSampler installation	Proceed to Connect the Pneumatics , Page 103.
you have been using your AutoSampler for Liquid Injection (and have therefore already connected the pneumatics)	Proceed to Remove Solvent from the Solvent Lines , Page 104.

CONNECT THE PNEUMATICS

The 8200 AutoSampler requires clean, chromatographically pure compressed (40 to 60 psig) air or nitrogen to operate using Ambient Headspace Injection. A chromatographic grade gas is recommended.

To connect the pneumatics, follow this procedure:

1. Connect the supply cylinder to the “barbed” fitting on the rear panel of the AutoSampler using the 1/8-inch OD polyurethane tubing from the Standard Accessory Kit. The barbed fitting is labeled GAS INLET 40-60 psig.



2. To make the connection at the GAS INLET, press the 1/8-inch OD tubing onto the barbed hose connector. Tug firmly on the tubing to make sure it is fully gripped on the fitting, referring to the diagram below.




WARNING: EYE HAZARD

Wear Proper Eye Protection. To prevent possible injury from a flailing gas line, shut off the gas supply valve before disconnecting the line.

NOTE: The pressure regulator has been preset at the factory to 20 psi.

1. Proceed to Page 105, **Install the 100- μ L Syringe Assembly.**



REMOVE SOLVENT FROM THE SOLVENT LINES

During the headspace sampling sequence of the 8200 AutoSampler, the solvent lines and syringe are flushed with air to minimize the carryover of sample between runs. Therefore, at least one of the solvent reservoirs must contain only air. If the unused solvent reservoir contains a solvent to which the detector to be used in the headspace analysis responds, make certain that you have removed all traces of this solvent, otherwise the solvent will interfere with the headspace analysis.

NOTE: You do not need to clear the solvent lines if the AutoSampler has not been used since it was shipped from the factory.

To remove solvent from the Teflon® solvent lines, proceed as follows:

1. Remove the solvent reservoir(s) that contain wash solvent from the AutoSampler.
2. Remove the solvent filter, nut and the blue ferrule from the Teflon® tubing. Carefully wipe any solvent from the outside of the tubing and let the tubing dry.
3. Replace the solvent reservoir with either a new solvent reservoir (P/N 03-908877-00) or the current solvent reservoir after it has been thoroughly cleaned and dried. You do not need to install a solvent filter. If a solvent filter is installed, replace that filter with a new solvent filter (P/N 28-211512-00) that has never been immersed in solvent.
4. Install a clean septum (P/N 00-997628-01) in the waste cup of the AutoSampler.

FLUSH THE SOLVENT LINES

To flush the solvent lines, proceed as follows:

1. Follow the steps given on **How to Inject a Single Sample** on Page 51, using the method you created in **Save Your Method** on Page 102. This procedure will cause the pure air or nitrogen to flush through the solvent lines even if you have not yet completed the remaining steps in this Ambient Headspace injection method setup.
 - Pure air or nitrogen is flushed through the solvent lines.
2. Repeat Step 1 at least three times to thoroughly flush the solvent lines.

INSTALL THE 100- μ L SYRINGE ASSEMBLY

Once you have correctly mounted the syringe in the syringe bracket, you will not need to make further adjustments.

NOTE: Do not install the 100- μ L syringe in the AutoSampler unless the solvent lines are free of solvent. If you need to clear the solvent lines, refer to **Remove Solvent from the Solvent Lines** on Page 104.


To install the syringe/needle assembly (P/N 03-918728-00) in the syringe carriage, proceed as follows:

1. Place both thumbs on the side locking tabs on either side of the syringe bracket. Push up on the bracket until the red tab appears. Open the front door of the syringe mount assembly.
2. Guide the syringe needle into and through the small hole in the needle guide. Press the syringe nut into the slot such that the top of the nut rests on the upper surface of the syringe carriage bracket. Press the syringe plunger button into the slot of the lead screw nut.
3. Close the front door of the syringe mount assembly and press the side locking tabs down while supporting the carriage. The red signal tab is no longer visible. With the locking tabs down, the syringe is locked into position so any movement of the syringe hex nut will be minimal.
4. Connect the solvent line to the syringe side arm. Slide the knurled nut back from the end of the tubing. Press the tubing firmly against the machined shoulder of the side arm port. Slide the nut up to the port and snugly thread it onto the side port. Tighten the nut with your fingers only. Do not use pliers.

ALIGN THE SYRINGE NEEDLE WITH THE INJECTOR PORT

To align the 100- μ L syringe needle with the injection port, proceed as follows:

1. Make sure that the AutoSampler power switch is OFF. The power switch is located on the rear of the tower.
2. Press down on the storage module latch release and swing the module out to access the injector assembly. With the storage module swung out and AutoSampler power switch OFF, a contact in the latch assembly is activated that toggles the carriage stop solenoid on or off each time the carriage position sensor is interrupted.



NOTE: To re-latch the Storage Module, the latch release must be pressed down.

3. Turn the GC power switch ON. Turn the AutoSampler power switch ON. Raise the syringe carriage from its rest position until you here a click. The stop solenoid toggles to permit you to lower the carriage by hand.

NOTE: As you align the syringe, take care to support the syringe carriage so that it does not drop and damage the syringe or needle assembly.

4. Loosen the screws holding the mounting plates in place so that you can move the AutoSampler on the GC cover. Lower the syringe carriage by hand until the needle tip is just above the GC injector needle port. Move the AutoSampler such that the syringe needle is precisely aligned with the injector needle port. Tighten the mounting plate screws securely.
5. Lower the syringe carriage two or three times so that the syringe needle punctures the injector septum. This ensures that the syringe and syringe needle are properly aligned with the injector port. Repeat this procedure each time you replace the injector septum.

NOTE: Once the mounting plate(s) have been aligned, they should not require further adjustment.

6. Raise the syringe carriage to its rest position (its upper limit of travel). Then, lower the syringe carriage a fraction of an inch so that it rests on the stop solenoid latch.
7. Close the storage module.

FINAL AMBIENT HEADSPACE SETUP CHECKLIST

Setup for Ambient Headspace Injection for your 8200 AutoSampler is now complete. Before operating the AutoSampler as described in Section 10, make sure you have

- Built the Ambient Headspace injection method
- Removed any solvent from the solvent lines
- Connected the pneumatics
- Installed the 100- μ L syringe
- Aligned the syringe needle with the injector port

Section 10

Ambient Headspace Injection Operation

This section gives additional information that you will need to operate the 8200 AutoSampler using the Ambient Headspace Injection Method. After familiarizing yourself with the information in this section, proceed directly to Section 5, **Performing Injections**.

NOTE: The Ambient Headspace Injection does not support liquid injection modes.

SAMPLING

LOADING THE CARROUSEL

The carousel rotates clockwise and the vial position numbers increase counterclockwise around the carousel. Observe this order when loading sample vials. Load sample vials according to the vial numbers shown on the carousel label.

The AutoSampler ignores vials without a septum. Seal all vials in the carousel that are to be sampled with a cap and septum.

100- μ L SYRINGE HEADSPACE INJECTION METHOD

The 48-vial carousel takes a standard 2-mL glass screw top vial with outside dimensions of 12 x 32 mm (12 x 35 mm with cap). The vial is sealed with a plastic screw cap with a 5-mm hole, and a white silicone rubber septum (8-mm diameter) faced with an inert film on one side. The inert side is the colored side of the septum. Place the colored side of the septum toward the sample solution.

- **Selecting Sample Vials.** A variety of sample vials are available from a number of sources. However, the vials you use with the AutoSampler must have the same external dimensions as the vials available from Varian (see **8200 AutoSampler Parts and Supplies**).
- **Installing Vial Septa.** The 2-mL vials (P/N 03-949835-00) are shipped complete with plastic screw caps and Teflon®-laminated silicone rubber septa. Always place the colored, chemically inert Teflon laminate face down (next to the contents of the vial).

Section 11

Automated Large Volume Injection Setup and Operation

INTRODUCTION

The Large Volume Injection (LVI) technique can improve GC and GC/MS detection limits by a factor of 10-50 times. This is done by transferring 10-50 times more analyte onto the GC column than with traditional injection, while eliminating the interference effects of excess solvent.

The LVI technique uses the 1079 Universal Capillary Injector and a 100 μ L liquid injection syringe. The sample is injected at a very slow rate while the injector temperature is set a few degrees below the solvent boiling point. The injector split valve is left open for a period of time to vent the solvent to the split vent. The split valve is then closed and the injector is rapidly heated to vaporize the solute material onto the GC column where the separation is made.

NOTE: Automated Large Volume Injection requires removal of the SPME agitation accessory. If this device is installed on your AutoSampler, refer to Section 8, Changing From SPME III to Liquid Sampling on Page 94.

COMPONENTS NECESSARY FOR LVI

1. 3800 GC
 - 8200
2. 1079 Universal Capillary Injector. Cryogenic cooling is recommended.
 - 2 mm insert packed with deactivated glass wool (P/N 03-925350-00), or
 - 3.4 mm insert packed with 10% OV-101 on Chromosorb W-HP (P/N 03-918956-00), or
 - 0.5 mm open insert (P/N 03-925331-00)
3. Star Chromatography Workstation with 8200 Control.
4. 100 μ L liquid syringe (P/N 03-925414-01).


CONDITIONS FOR LVI WITH VARIOUS SOLVENTS

The following tables give all of the conditions for making a 10-50 μL injection of a sample dissolved in a low boiling solvent (in this case MeCl_2 , bp 42°C) and a higher boiling, higher viscosity solvent (in this case EtOAc , bp 77°C). The GC column, column ramp rate, and final temperature will have to be chosen based on the compounds you are separating. The conditions given here are for anthracene. You can use these examples to develop the conditions for your LVI application. A detailed description of the LVI procedure has been included so that you can understand the purpose of the selected times, flow rates and temperatures.

The initial injector temperature is set a few degrees below the boiling point of the sample solvent. This prevents the solvent from vaporizing too quickly which might cause a loss of sample through the split vent. The injector split vent is set to 100 times the column flow so as to efficiently remove the solvent from the injected sample. If the split vent flowrate is set too low, excessive solvent will be injected onto the column resulting in peak shape distortion.

After 1.5 minutes the split vent is closed and the temperature of the injector is rapidly raised to 325°C . The sample is then vaporized onto the column which is held at the same temperature as the initial injector temperature. At 3.0 minutes the split vent is opened to vent any residual sample or solvent, cleaning the insert for the next injection. The injector is kept at 325°C until 3-5 minutes before the end of the method, at which time the injector is programmed to its initial temperature. If cryo cooling of the injector is used, you will have to adjust the time at which you program the injector to return to its initial temperature. Make it a target to have the injector and column oven reach their initial temperature at about the same time. This is a way to shorten the time between runs and minimize use of coolant.

Between 1.5 and 3.0 minutes the sample (less most of the solvent) is moved to the head of the column by closing the split vent and heating the injector. At the same time that the split vent is opened at 3.0 min the column temperature ramp is started. The column ramp is started to separate the components. See Diagram 1 for a graphic description of the LVI procedure.



The length of time the split vent is closed depends on the injection volume. The smaller the injected volume the less time will be required to vent the solvent. For example, with a 10 μL injection the split vent can be closed after 1 min instead of 1.5 min. You would also need to adjust the time the injector temperature ramp is started to 1.0 min as it is synchronized with the closing of the split vent. Usually it is necessary to keep the split vent closed for 45-90 seconds depending on the amount of sample injected. The larger the sample the longer the time. The 90 second time was used successfully for an injection of 50 ng of anthracene.

The 8200 AutoSampler must be used in the User Defined mode so that four critical parameters can be adjusted. These parameters are solvent plug size, uptake speed, injection rate, and pause time. Solvent flush sampling is selected with both upper and lower air gaps chosen. The solvent plug size is set to 5 μL . Setting it larger can cause the upper air gap bubble to move to the plunger tip mixing the sample with the solvent plug. If this occurs the sample will not be completely flushed out of the needle. The uptake speed and injection rate are set to their lowest values. The pause time is set to its maximum (10 seconds) to allow sufficient time for the syringe to completely fill after sample withdrawal. If the sample is withdrawn too quickly cavitation will occur and more than 10 seconds (the maximum pause time) will be required for sample to fill in the cavity (bubble). Also, if injection rate is too fast the sample components might be lost by being carried out the split vent with the high concentration of solvent vapor generated in the insert.

Suggested Injector Inserts for LVI

- 2 mm packed with deactivated glass wool (03-925350-00). This insert is used with relatively clean samples. It is a general purpose insert shipped with the injector.
- 3.4 mm packed with 10% OV-101 on Chromosorb W-HP 80/100 (03-918956-00). This insert is useful for dirty samples and provides added retention for your sample components.
- 0.5 mm open insert (03-925331-00). This insert is recommended for trace levels of polar compounds where maximum inertness is desired. The recovery of sample with this insert can be significantly less than for the two above inserts.

The maximum recommended temperature for the above three deactivated inserts is 325°C. If deactivation is not an issue, higher temperatures can be used. The 10% OV-101 packing is rated to 350°C.

Additional Hints and Procedures

1. **100 μ L syringe installation.** The procedure for installing the 100 μ L syringe is the same as the 10 μ L syringe. Be sure to adjust the syringe bracket position so that there is no gap between the end of the syringe needle nut and the bracket it rests on.
2. **Plugged needle.** If your needle becomes plugged during use it will most likely cause a stall of the syringe plunger stepper motor causing the syringe to stop with the needle inserted in the waste cup. If you find your system has stopped and the needle is in the waste cup, you have a plugged needle.
3. **Replacement of needle.** Unscrew the nut holding the needle in the syringe. Pull the needle out of the syringe body. Carefully insert the new needle and tighten the nut 1/8-1/6 turn after the nut first starts to tighten. The needles for the 10 and 100 μ L syringes are different. The needle for the 100 μ L syringe uses a black Kalrez® seal without a spring while the 10 μ L syringe needle uses a white Teflon® seal with a spring.
4. **The plunger seal** at the top of the syringe is adjusted at the factory and will generally not need adjustment.
5. **Syringe leak test.** To test the syringe plunger for leaks first remove the syringe from the 8200. Pour a few mL of solvent (whatever you are currently using) in a 25-50 mL beaker. Place the needle in the solvent and withdraw some solvent into the syringe. Withdraw the needle from the solvent and push the syringe plunger at a moderate rate. If the plunger is leaking solvent

will squirt out of the side arm. If there is no leak, the solvent will squirt out of the needle.

6. **Suggestion for Data Handling.** The solvent front for LVI is larger than with regular splitless injection. Setting Inhibit Integrate (II) in the data system timed events from 0-7 min prevents peak processing of the solvent peaks.
7. **Solvent purity.** Since you are injecting large volumes of solvent, you are now concentrating the impurities of the solvent. With LVI it is important to use high purity solvents for samples and wash solvents to minimize extra peaks in your chromatographic run. Run a blank for your solvent so that you know what peaks it is adding to your chromatogram.

Table 11-1: 8200 Star Workstation AutoSampler Conditions for LVI

	Star Chromatography Workstation Setting
Sample Mode	User defined
Syringe	100 μ L
End time	1.00 min
Solvent	A or B (same as sample solvent)
Solvent flush sampling	Yes
Syringe wash	30 sec
Air dry after wash	No
Solvent plug size	5 μ L
Vial needle depth	90%
Uptake speed	10 μ L/sec
upper air gap	Yes
lower air gap	Yes
Pause time	10 sec
Hot needle time	0.00 min
Injection rate	2.0 μ L/sec
Needle residence time	0.05 min

Table 11-2: 3800 GC Conditions for LVI

	For EtoAc (bp 77°C)	For MeCl ₂ (bp 42°C)
1079 INJECTOR		
Initial temperature	75°C	40°C
Initial hold time	1.50 min	1.50 min
Final temperature	325°C	325°C
Rate	200°C/min	200°C/min
Hold time	10.90 min	8.90 min
Coolant to injector on	Yes	Yes
GC COLUMN OVEN	—	—
Initial temperature	75°C	40°C
Initial hold time	3.00 min	3.00 min
Thermal stabilization time	2.00 min	2.00 min
Final temperature	270°C	270°C
Rate	15°C/min	25°C/min
Hold time	2.00 min	3.80 min
SPLIT STATE (EFC)	—	—
Initial Split state	On	On
1.50 min	Off	Off
3.00 min	On	On
COLUMN/INJECTOR FLOW RATES/MISCELLANEOUS		
Column	30 meter 250 μ ID 0.25 μ film DB-5	30 meter 250 μ ID 0.25 μ film DB-5
Column flowrate	1.5 mL/min at 75°C (16 psi He)	~1.8 mL/min at 40°C (16 psi He)
Injector septum purge	5 mL/min	5 mL/min
Split vent flow rate	150 mL/min	150 mL/min
Test sample	Anthracene (2 ng/ μ L in EtOAc)	Anthracene (2 ng/ μ L in MeCl ₂)
Wash Solvent	EtoAc	MeCl ₂


Section 12

Maintenance/Troubleshooting

Routine and preventative maintenance procedures are described in this section. When followed, they ensure optimum performance of the 8200 AutoSampler. Varian does not recommend a specific schedule for these maintenance procedures. Rather, you should determine your own regular maintenance schedule for the AutoSampler that reflects the frequency and duration with which the AutoSampler is used. Pages 117 through 121 describe specific maintenance operations which may be performed by the user. Refer to the following diagram.

NOTE: Lubricate or clean only those parts referred to in this Maintenance/Troubleshooting section.

- A Syringe Plunger Lead Screw
- B Belt Tensioner/Pulley
- C Syringe Carriage Stop Latch Solenoid
- D To Position #4 on J8
- E Syringe Carriage Drive Belt
- F To Position #7 on J8
- G Syringe Carriage Drive Motor
- H Syringe Carriage



SYRINGE CARRIAGE GUIDE RODS

The guide rods and syringe carriage bushings do not require lubrication. However, you should periodically wipe the guide rods with a clean, lint free cloth dampened with isopropanol. How often the syringe carriage guide rods need cleaning is determined by the cleanliness of the environment in which the 8200 AutoSampler is used.

SYRINGE CARRIAGE DRIVE BELT


The drive belt is adjusted at the factory and does not usually require adjustment. If you need to adjust or replace the drive belt, proceed as follows:




WARNING: **SHOCK HAZARD**

**Dangerous Voltages Exposed When High Voltage Cover is removed.
Unplug Power Cord.**

1. Turn off the power to the AutoSampler and disconnect the power cord. The power switch is located on the rear of the tower.
2. Remove the module cover by removing four screws: the two screws on the top of the cover and the two screws at the lower rear panel.



3. Remove the four screws at the module connecting brackets.
 - Disconnect the flat cable connector from J9 on the Controller PC Board and position #4 connector of J8. See the following diagram.
 - Separate the injection tower module and the pneumatics module.
4. Locate the belt tensioner bracket and mounting screws. Loosen the bracket mounting screws.
 - Adjust the bracket to increase/decrease the tension on the belt. There should be just enough tension to prevent slippage of timing belt teeth on the motor drive timing belt sprocket.
 - Tighten the tensioner bracket screws.
5. Reconnect the position #4 connector of J8, and connect the flat cable connector to J9.
 - Replace the tie wrap on the connectors if removed. Reassemble the modules and tighten all screws.
6. Reconnect the power cord.



SOLVENT COMPATIBILITY

The materials in the 8200 AutoSampler that contact solvents are passivated 316 stainless steel, inert polymers and fluoropolymers. Any solvents used with the 8200 AutoSampler must be compatible with these materials.

Avoid the following solvents:

- Halide- (F, Cl, BR, I) containing buffers or acids at concentrations exceeding 2M or having pHs of less than 2.
- Solvents that can form HCl, such as carbon tetrachloride or chloroform.
- Freon 11, TF, or FC-72.

SYRINGE AND SYRINGE NEEDLE MAINTENANCE

SYRINGE CARE

Acetone is recommended as a cleaning solvent when cleaning the syringe. Do not immerse the syringe in the solvent. The adhesive used to bond the metal flange and front needle fittings to the glass barrel will be damaged by soaking in a solvent. Clean only the outside surfaces of the syringe by wiping with a cleaning material.

Thoroughly rinse with acetone and dry before storing the syringe.

Varian syringes have been designed to withstand temperatures to 120°C. Rapid changes in temperature, however, may cause the glass barrel to split.

NEEDLE CARE AND REPLACEMENT

Needles are checked for correct flow before shipping. If your needle does become blocked during usage, remove it from the barrel and gently push solvent through it in the reverse direction.

Check the needle tip for burrs that can damage septa material, causing septa material to deposit in the chromatographic system.

If you need to replace the syringe needle, you will need

- Replacement needle (P/N 03-918987-00)
1. Remove the syringe/needle assembly from the syringe mounting bracket.
 2. Unscrew the syringe nut and remove the needle assembly.
 3. Remove the syringe nut and spring from the old needle; fit them to the new needle as indicated in the diagram.
 4. Introduce the back of the needle into the syringe barrel. Push in on and tighten the syringe nut with your fingers.
 5. Snap the syringe/needle assembly into the syringe mounting bracket. Check that the needle is properly positioned over the injector port.

REPLACING THE MAIN POWER FUSE


To perform this operation, you will need a narrow blade screwdriver.

To replace the main power fuse, proceed as follows:


1. Open the fuse block/cover assembly with a narrow blade screwdriver.
2. Snap out the defective fuse and replace with a fuse of the same voltage and amperage.

NOTE: One 0.5 Amp SB fuse will be exposed for 100/120V. Two 0.2 Amp T-Type fuses will be exposed for 220/240V.

Fusing Arrangement (220/240V)



Fusing Arrangement (100/120V)



INSTRUMENT TROUBLESHOOTING

Refer to Table 11-1 for descriptions of possible chromatographic and hardware problems, their causes, and the actions you should take to correct them.

NOTE: Before troubleshooting a problem with your 8200 AutoSampler, make one or more manual injections first to establish that the problem is in fact a problem with the AutoSampler. Often, the source of the problem is misdiagnosed as a problem with the AutoSampler when the problem is actually with the gas chromatograph (leaks, poor column/injector connections, adsorption within the injector).

Table 12-1 Troubleshooting

Observed Problem	Possible Cause	Corrective Action
Split peaks	Injection speed too slow for split injection	Use a fast injection speed (5 to 10 $\mu\text{L}/\text{sec}$) for split injection to obtain best peak shape, especially for early eluting peaks.
	Injection speed too fast for solvent/column combination	Narrow bore capillary columns easily produce split peaks due to solvent flooding. Reduce injection speed and raise initial column temperature to vaporize solvent gradually during cold on-column or splitless injection.
	Polar sample solvent beads up on non-polar or semi-polar liquid phase	Use a non-polar solvent plug with polar sample solvents to act as co-solvent.
	Sample size too large	Can also cause solvent flooding with narrow bore columns. Reduce sample size if possible; reduce injection speed. Reducing solvent plug size may also help, but may lead to increased solute discrimination.
	Component concentration too high	Exceeding the capacity of the column liquid phase (overloading the column) may cause split or deformed peaks. Reduce the sample concentration by diluting the sample.
Poor peak area reproducibility (High Area RSDs)	Syringe leaking at needle seal (solvent appears around syringe nut)	Check that the syringe nut is tight; if leaks continue, replace needle and Kalrez® seal. When replacing needle, check syringe barrel surface for cracks or chips. If you notice cracks or large chips, replace the syringe.
	Syringe leaking at plunger seal (solvent appears on upper body of syringe)	Leaks at plunger seal lead to solvent flowing over outside of syringe barrel during wash cycle. Tighten plunger guide to restore seal. Do not overtighten; increased friction causes stepper motor to stall.
	Syringe plunger tip Teflon® seal worn	Replace syringe. A worn Teflon plunger tip seal allows air or solvent to flow past the plunger tip, causing inconsistent sample sizes to be withdrawn from vial.
	Viscous sample injected with wrong mode	Use Viscous Sample Mode, instead of Standard Mode.

Observed Problem	Possible Cause	Corrective Action
	Volatile sample injected with wrong mode	Use Volatile Sample Mode instead of Standard Sample Mode when injecting solvents such as methylene chloride, carbon disulfide, or pentane.
	Partially plugged needle	Replace needle if you cannot clear the plug.
	Inlet problems	Check for leaks, flow controller failures, pressure fluctuations, etc. See Troubleshooting section of your GC Operator's Manual.
	Waste arm plugged	Check for free flow of wash solvent through syringe and out waste arm.
	Solvent reservoir empty	Refill reservoir.
	Too many injections per vial	With long run times and volatile sample solvents, concentration of components in sample increases after sample vial septum is punctured. When using Neat Mode, a slight vacuum forms within vial as sample is withdrawn.
	Insufficient solvent plug volume when using solvent flush injection	Solvent plug volume must be greater than needle volume; i.e., 1.0 μ L or more.
Bent needle	Tray misaligned	If syringe is installed correctly in syringe holder (vertical, with no binding), and needle hits edge of vial, call Varian Customer Support. Do <i>not</i> rotate the carousel to achieve proper alignment!
	Wrong needle insertion depth	With certain disposable vials and some microvials, a needle depth default of 90% allows needle to hit bottom of vial. Change needle depth to lower value (e.g., 80%).
	Wrong vial used	Use only 1.5 mL vials and microvials certified acceptable by Varian. In particular, some microvial inserts have very small openings, allowing the AutoSampler needle to hit the glass edge.
	Syringe misaligned	Adjust the syringe mount, following the procedure in described in this Operator's Manual.
	Needle hits edge of septum nut when injecting	If syringe is mounted correctly in syringe holder (vertical, with no binding, aligned with sample vials), then loosen the AutoSampler mounting plate screw and align the AutoSampler with injector, as described in this Operator's Manual.

Observed Problem	Possible Cause	Corrective Action
	Too much resistance from new septum	After a new septum is installed, lower the AutoSampler carriage by hand until the needle pierces the new septum. This procedure also checks the alignment of the needle with the injection port. Do not overtighten the injector nut, as this increases the force required for the needle to penetrate the septum.
	Obstruction within injector	Check for obstructions or for incorrect column installation in injector. In particular, make sure that small diameter columns are not being used with large diameter inserts in the 1079. The wrong insert allows columns to be inserted too far into the injector.
Sample component carryover ("ghost peaks")	Wash time too short	Use longer wash time (default of 20 seconds).
	Wrong wash solvent	Use a wash solvent that is miscible with sample solvent. If the sample is composed of constituents with a wide range of polarities, especially in high concentrations (e.g., perfumes), use two wash solvents of different polarities.
	Defective syringe and/or needle	Voids or gaps in syringe seals or dead volume in tip of defective needle leads to carryover; replace syringe and/or needle.
	Inlet problem	Inlet may be too cold, poor seals may produce poorly swept areas, or otherwise contribute to carryover. Consult Troubleshooting Section of your GC Operator's Manual.
	Plugged wash bottle filter	Replace filter. Verify that wash solvent flows freely through syringe during wash cycle.
	Valves plugged due to swollen Viton® seals caused by using ketones as wash solvent	If necessary to use ketones as wash solvents, replace standard valves with optional valves built with Kalrez seals (P/N 03-918533-90).
	Waste arm plugged	Clear waste arm; ensure free flow of solvent through syringe during wash cycle.
	Wash solvent reservoir empty	Refill reservoir.

Observed Problem	Possible Cause	Corrective Action
	Contaminated gas supply to the AutoSampler	Not actually sample component carryover, but extra peaks due to contaminants dissolved in wash and backup solvent. Check that the gas supply to the AutoSampler is chromatographically pure. At minimum, direct gas supply through a good hydrocarbon/organics filter.
	Contaminated wash solvent	Use only chromatographic grade wash and backup solvents.
	Septum contamination	Condition injector septa in GC oven prior to use. For ultratrace analysis, condition waste arm septa also.
Wash solvent carryover	Wrong sampling mode chosen	If wash solvent carryover is not desired, use the Neat Sampling Mode.
	Wash solvent too viscous to be removed by Neat Sampling Mode	Use a less viscous wash solvent; or use the original wash solvent followed by a less viscous solvent; or use air dry in User Defined Mode to remove most solvent before sampling.
	Wash bottle filter and lines not purged properly when changing to new solvent	When changing solvent in a wash bottle, use a new, dry filter frit. The porous frits absorb a large quantity of solvent, which is difficult to remove completely by purging. After changing solvents, set wash time to 180 sec and repeat several wash cycles.
No peaks, or small peaks	Plugged syringe needle	If you cannot clear the syringe needle, replace needle.
	Bad syringe needle seal	A poor needle seal allows air to be sucked into syringe instead of sample. Replace needle and seal; ensure syringe nut is tight.
	Syringe needle not inserted into sample	Increase value of Vial Needle Depth so that syringe needle dips into the sample. Carefully adjust the Vial Needle Depth when sampling volumes of 20 μ L or less.
	Teflon® plunger seal worn or defective	Sample may be withdrawn into syringe, but leaks past defective Teflon seal when needle is inserted into pressurized injector. Replace syringe.
	GC problems	Wrong splitter setting, leaks, wrong detector settings, broken column, etc. Consult Troubleshooting Section of your GC Operator's Manual.

Observed Problem	Possible Cause	Corrective Action
Short Syringe Lifetime	High friction due to use of air dry	Use air dry only if absolutely necessary.
	Overuse of Neat Sampling Mode	Use the Neat Sampling Mode only when you need to minimize the wash solvent carryover. Syringe wear is increased when using this mode.
Syringe needle pulled vial out of tray	Vial retainer arm misaligned	Realign retainer arm or call Varian Customer Support.
	Unapproved vial in carousel	Use only vials approved for use in the 8200 AutoSampler.
Vial not found or skipped	Cap of vial removed or missing	Infrared sensor beam must be blocked to detect vial. If vial must be sampled without cap, place opaque tape over the bottom of the vial.
	Defective infrared Sensor	Call Varian Support.
Solvent or sample peak tailing	Injection speed too fast for splitless or large bore (0.53 mm) vaporizing injection	Use slow injection speed (1-2 $\mu\text{L}/\text{sec}$), especially with large sample volumes.
	Inlet or other GC problems	Poor internal seals, plugged septum purge, too low a vent flow, etc., contributes to tailing. Consult Troubleshooting Section of your GC Operator's Manual.
Poor retention time reproducibility	Rarely due to AutoSampler malfunction	Often an inlet system problem. Change injector septum first. If problem persists, consult Troubleshooting Section of your GC Operator's Manual and check entire system for leaks. If you suspect an AutoSampler timing problem, call Varian Customer Support.
Sample discrimination (early or late eluting peaks reduced in size relative to other components)	8200 AutoSampler not set to Solvent Flush Mode	Standard Mode, which uses solvent flush injection, should be used to minimize sample discrimination. Do not change the default parameters.
	Injection speed too high	Too high an injection speed into a hot injector results in flashback and the loss of low boiling components. Reduce the AutoSampler injection speed.
	Wrong choice of backup solvent	Low boiling solvents vaporize too rapidly from the syringe needle, leaving higher boiling components behind, especially with slow injection speeds.

Observed Problem	Possible Cause	Corrective Action
	Wrong injector temperature	Too high an injector temperature can cause low-boiling solutes to be lost. Too low an injector temperature can cause high boiling components to be lost in a vaporizing injector. Ideal results are obtained with non-vaporizing injection (1079), followed by temperature programming of the injector.

Section 13

8200 Parts and Supplies

Description	Part Number
10 µL Liquid Injection Syringe	03-918986-00
10 µL Syringe Needle	03-918987-00
Interface Cable 3800 to 8200	03-925302-01
SPME Test Sample	03-918967-00
12-Vial Carrousel	03-918870-00
Labels for 12-vial Carrousel	03-925270-00
Labels for 48-vial Carrousel	03-925271-00
10 mL Cap: Vials (pk/36)	03-918873-01
10 mL Cap with Septum (pk/36)	03-918875-01
O-ring for 10 mL Vial (pk/36)	03-918874-01
10 mL capacity Vials (pk of 288)	03-918873-02
10 mL cap with Septum (pk of 288)	03-918875-02
O-ring for 10 mL Vial (pk of 288)	03-918874-02

Large Volume Injection	
100 µL Syringe with Needle	03-925414-01
100 µL Syringe Needle	03-908824-00
2 mm 1079 insert packed with deactivated glass wool	03-925350-00
3.4 mm 1079 insert packed with 10% OV-101 on Chromosorb W-HP 80/100	03-918956-00
0.5 mm 1079 open insert	03-925331-00

Ambient Headspace	
8200 Syringe, 100 µL w/sidearm	03-918728-00
100 µL Syringe Needle	03-908824-00

SPME FIBERS (3 FIBERS PER PACKAGE)

Polydimethylsiloxane

100 μm coating for volatiles (280°C inlet, red hub)

Description	Part Number
For manual sampling	03-918963-01
For 8100/8200 AutoSamplers	03-918963-02

30 μm coating for semivolatiles and pesticides (280°C inlet, yellow hub)

For manual sampling	03-918963-09
For 8100/8200 AutoSamplers	03-918963-10

7 μm coating for mid to nonpolar semivolatiles (320°C inlet, green hub)

For manual sampling	03-918963-04
For 8100/8200 AutoSamplers	03-918963-03

Polyacrylate

85 μm coating for semivolatiles (310°C inlet, white hub)

For manual sampling	03-918963-05
For 8100/8200 AutoSamplers	03-918963-06

Carbowax/Divinylbenzene

65 μm coating for alcohols and polar compounds (265°C inlet, orange hub)

For manual sampling	03-918963-11
For 8100/8200 AutoSamplers	03-918963-12

Assortment Pack

One fiber of each type (100 μm PDMS, 7 μm PDMS, 85 μm Polyacrylate)

For manual sampling	03-918963-07
For 8100/8200 AutoSamplers	03-918963-08

SPME Holders

For manual sampling	03-918964-01
For 8100/8200 AutoSamplers	03-918964-02

Section 14

Application Notes

To obtain copies of the following Application Notes, or for other information about Varian Products, contact your local Varian Sales Representative.

- Method Development Tips for the Automated SPME System (Varian GC Advantage Note 11) (Included in this section)
- Determination of Halogenated Pesticides with Automated SPME and Agitation (Varian GC Advantage Note 17)
- A Demonstration of Automated Solvent Flush Injection to Eliminate Sample Discrimination in the Determination of pesticides by Splitless Capillary Injection (Varian Application Note 28)
- Automatic Injection of Difficult Samples: Fine-tuning the Conditions for Viscous Samples (Varian Application Note 42)
- Automatic Injection of Difficult Samples: Fine-tuning the Conditions for Volatile Samples (Varian Application Note 43)
- Determination of Volatile Organic solvents in Water by Headspace Sampling with the 8200 AutoSampler (Varian Application Note 46)
- Profiling Flavors in Alcoholic and Non-Alcoholic Beverages with Automated Solid Phase Microextraction (SPME Application Note 1)
- Determination of Residual Solvents in Pharmaceuticals with Automated Solid phase Microextraction (SPME Application Note 2)
- Determination of a Wide Range of Organic Impurities in Water with Automated Solid Phase Microextraction (SPME Application Note 3)
- Flavor Analysis of a Fruit Beverage with Automated Solid Phase Microextraction (SPME Application Note 4)
- Analysis of Thermanol in Process Water Using Solid Phase Microextraction (SPME Application Note 5)
- Characterization of Flavor Components in Wines with Solid Phase Microextraction (SPME), GC and GC/MS (SPME Application Note 6)

- Determination of Residual Solvents and Monomers in Polymers with Solid Phase Microextraction (SPME) and GC/MS (SPME Application Note 7)
- Determination of Trace Methanol in Caustic Industrial Product with Automated Solid Phase Microextraction (SPME Application Note 8)
- Blood Alcohol Determination with Automated Solid Phase Microextraction (SPME Application Note 9)
- Rapid Analysis of BTEX and TPH in Water using Solid Phase Microextraction (SPME Application Note 10)
- The Advantages of SPME Agitation (SPME Application Note 11)

Method Development Tips for the Automated SPME System

GC Advantage Note 11

Zelda Penton, Varian Chromatography Systems

Introduction

Automated solid phase microextraction (SPME) can yield ppb detection limits for organic compounds in water or solids. Linearity is excellent, and relative standard deviations are often better than 3%. However, the technique is new and many questions have been raised by chemists who are excited by the possibilities and would like to utilize SPME in their analytical procedures. The questions fall mostly in the category of SPME hardware and sampling parameters. These questions will be addressed here. Experience has shown that sample preparation is the key to good results with SPME; therefore, techniques for working with volatile samples will also be discussed.

For recommendations on injectors and injector inserts, see the SPME operator's manual. Some guidelines to help the user get started are given below. These suggestions are discussed in greater detail in the ensuing sections.

Guidelines For Getting Started

1. A new fiber must be conditioned by desorbing for a minimum of five minutes in an injector that is at least 10°C hotter than the temperature to be used during the analysis. The GC column should then be temperature programmed; this procedure should be repeated until there are no extraneous peaks.
2. The fiber and injector septum should be changed after approximately 100 runs.
3. For water soluble analytes, sensitivity can be enhanced by saturating the sample with salt (usually sodium chloride or sodium sulfate).
4. Standards and samples are normally prepared and diluted in **storage containers** and then transferred to AutoSampler **vials** for analysis. Prepare samples and standards carefully so that volatiles are not lost:
 - a. After preparation, water samples containing volatiles should completely fill the storage container without any headspace.
 - b. Store samples in the refrigerator. Chill the AutoSampler vials before adding the sample.
 - c. Transfer samples to the AutoSampler vials with a pipette of sufficient capacity to deliver the entire sample in one step. The outer diameter of the pipette should be small enough to allow the pipette to **easily fit** into the AutoSampler vial.
5. For headspace sampling, use up to 0.8 mL sample; for liquid sampling use 1.2 mL. Do not fill the AutoSampler vial to the top!
6. A reasonable adsorption time is 15 minutes with 2 minutes desorption; however, these conditions should be optimized for each analysis. It is not necessary to achieve equilibrium if the total analysis time will be prolonged. For many samples, RSD's under 5% can be obtained prior to reaching equilibrium.

Hardware And Software

What fibers are available?

The 100- μ m polydimethylsiloxane (PDMS) fiber should be used for the majority of applications but two other fiber types are available. The features of currently available fibers are summarized below.

PHASE	ADVANTAGES	DISADVANTAGES
100- μ m PDMS	High sample capacity. Suitable for a wide variety of applications- volatile compounds such as dichloromethane to semi-volatiles (some pesticides).	Longer adsorption time required to reach equilibrium. Possibility of carryover for some semivolatiles.
7- μ m PDMS	Bonded fiber which allows a higher desorption temperature than the 100- μ m PDMS (Maximum operating temperature 320°C vs. 260°C for the 100- μ m fiber). Faster desorption-therefore initial column temperature during the GC run can be higher. This shortens the total run time. Thinner phase means shorter time to reach equilibrium.	Low sample capacity. Not suitable for volatiles.
Polyacrylate	Greater affinity for polar compounds than above fibers.	This phase is more of a solid than the PDMS phases; therefore diffusion rates are slow and equilibration times are long.

All of the data presented in the figures and tables was obtained with 100- μ m PDMS fibers.

How long does a fiber last?

The fiber life will vary with experimental conditions but typically, there is no evident deterioration in chromatography up to 100 runs when desorbing into an injector heated to 220°C. This is true even when immersing the fiber into water that is saturated with salt and is at pH 2. One sign of an aging fiber is deterioration of precision. This might also be due to the aging of the septum. It is best to change the septum when changing the fiber.

Is it necessary to use a septumless injector?

Users sometimes express concerns that the sheath on the fiber rod might core the injector septum. Figure 1 shows no change of retention time after 46 runs, indicating that the septum is intact. Again, the practice at Varian is to change the septum and the fiber at the same time after about 100 runs according to the following procedure:

1. A new septum is installed.
2. The AutoSampler carriage is lowered manually so that the fiber sheath pierces the septum several times.
3. The new septum is removed and inspected. If loose particles of septum material are present, they are removed.
4. The septum is reinstalled.





Figure 1. Demonstrating the integrity of the injector septum after 46 SPME injections. Note that the scale on the Y axis is 0.04 minutes. A Thermogreen™ LB-2 (Supelco) septum was used in this study.

How much sample should be added to the AutoSampler vials?

The AutoSampler vial should not be filled to the top. It was observed when sampling volatiles, that equilibrium was attained faster when headspace was present, even when liquid was being sampled. Furthermore, immersion of the metal fiber-support rod in the liquid sample may result in the adsorption and/or breakdown of analytes.



A final reason not to fill the vial, is the possibility of carryover if liquid sample enters the fiber sheath. A sample volume of 1.2 mL is recommended for liquid sampling and the maximum volume should be 0.8 mL when sampling the headspace.

What are the recommended adsorption and desorption times?

Adsorption time varies inversely with the volatility of the analyte and also depends upon the relative volumes of the phases in the vial. Satisfactory precision can often be obtained without achieving equilibrium. This is convenient if the GC cycle time is relatively short and prolonged sampling times would greatly lengthen the total analysis time. A reasonable point for sampling is fifteen minutes, but adsorption times may be longer if the GC cycle time permits.

At least 2 minutes is recommended to desorb all traces of the analyte to minimize carryover. The injector temperature should normally be at least 200°C but should not be higher than the temperature limit of the analytical column or the SPME fiber. If carryover is present, a longer desorption time and/or higher injector temperature should be used.

Is cryofocusing necessary?

Injector: After initial studies, it was concluded that injecting into a hot injector gives the best results, even when the sample contains very volatile analytes such as vinyl chloride.

Column: With a 0.53 mm column coated with a 3-micron moderately polar phase, cryogenic focusing is required in the column oven to elute vinyl chloride with a good peak shape. For dichloromethane (BP 40°C) and compounds with higher boiling points, cryogenics are unnecessary. SPME sampling does not require special GC conditions.

Is the fiber easily saturated and do compounds tend to be displaced in mixtures?

This depends on many factors including sample size, the affinity of the fiber for the components of a particular sample, and the fiber thickness. It was found in recent studies that there was a linear response to benzene at concentration ranges greater than 300 ppm. In another experiment, benzene gave the same response whether it was the only organic compound in water or in a test mixture containing several other compounds (Table 1). Nevertheless, it is important in developing and validating a method, to do recovery studies. Displacement effects may be minimized by diluting the sample to a lower concentration. With 2-mL AutoSampler vials, the probability of fiber saturation is reduced.

conc. (ppm)	MeCl ₂	CHCl ₃	Benzene	TCE	Dioxane	Toluene	Xylene	TMB*
1	4281	3567	175183	51579	1203	308884	427923	637614
2	8471	6322	337571	101351	2165	617885	872322	1251551
4	17894	12478	704235	207027	4192	1233313	1713901	2331268
1			154458					
2			328038					
4			674635					

Table 1. FID area counts for benzene alone in water at (bottom 3 rows) and in a mixture with several other organic compounds in water (top 3 rows). All of the compounds were at the concentration shown at left except chloroform which was at half the concentration shown. These data were obtained by SPME fiber sampling of the liquid phase (1.2 mL).

*1,2,4-trimethylbenzene

Is it better to sample the liquid or the headspace?

Theoretically, the response should be the same if the volumes sampled are similar. This has been found to be the case for the compounds listed in Table 1. Practically, for compounds of very low volatility, the adsorption time from the headspace would be long and liquid sampling would be preferable.

Does the fiber pick up contaminants in the air that will interfere with the analysis?

A new fiber needs to be desorbed for 15-20 minutes. When a new 100-micron fiber is desorbed for 3 minutes, followed by a GC run, typically, after six GC runs, the blank is clean. For a fiber that has been used, the first run each day should be a blank. Ghost peaks occasionally appear from AutoSampler vial septa. If this occurs, the user should try a different brand of septum and/or bake the septa before use.

In order to minimize the presence of extraneous peaks, the SPME software parameters should be set so that the GC is ready for the sample to be injected immediately after adsorption. See the SPME manual for a detailed explanation.

Sample Preparation

Experience has shown that poor precision and accuracy often result from improper sample handling.

- It is a common practice to prepare a stock solution of analytes in methanol and then add a small aliquot to water. This method is acceptable with SPME and results are the same as adding organics directly to water if the total level of methanol is less than 1%.
- Saturation of the standards and samples with sodium chloride or sodium sulfate is useful in two situations:

The analytes are polar and soluble in water.

The samples contain salts and it is desired to minimize matrix differences.

- When determining acidic compounds such as phenols, lower the pH; for basic compounds such as amines, raising the pH will enhance sensitivity.
- It is important to keep the concentrations of all of the components in an aqueous solution sufficiently low so that they remain dissolved in water.


For volatile compounds, additional guidelines should be followed:

- When preparing standards of volatiles in water, the liquid should fill the entire storage container without any headspace.
- Losses can occur when diluting the high level standard in preparation for a linearity study. To minimize errors, fill the containers that are to contain diluted standards with cold water at the correct volume for the dilution, quickly pour the concentrated standard into the containers, cap, mix and refrigerate. For example, when diluting to 1/2 and 1/4 the concentration of the highest standard, take 40 mL vials (44 mL when filled to the top), add 22 mL and 33 mL of cold water; then pour in the concentrated standard.

- If salt is added or the pH is adjusted, great care should be taken to minimize losses. For example, if the samples and standards are to be diluted in water, the salt can be added to the water *before* the dilution is made.

Chill the AutoSampler vials before adding the samples. Remove the standards and samples from the refrigerator, uncap them and quickly transfer aliquots to the AutoSampler vials, using a pipette that easily fits into the neck of the vial.

Ideally, the sample will be transferred with a disposable 1-mL serological pipette that is calibrated to deliver up to 1.2 mL and is graduated in 0.1 mL increments.



- Cap the AutoSampler vials quickly. If solids were placed in the vial, prior to adding the liquid, mixing with a vortex mixer will assure a homogeneous sample. Aqueous standards and samples remaining in the storage containers that were used to fill the AutoSampler vials, should not be used again.

Just before analysis:

One or two blanks should be run. The AutoSampler vials containing the samples should be allowed to reach room temperature before sampling.

Index

.mth extension45, 52, 72, 102
100 μ L Headspace Syringe 99
100 μ L Syringe 101
12-Vial Carrousel.....5, 7, 8, 69, 78, 84, 129
48-Vial Carrousel.....5, 7, 8, 50, 83, 109, 129
8200 AutoSampler, Features..... 5

A

Absorb Time 71
Activate a Method..... 51, 60
Activator Strip 91, 97
Active SampleList.....56, 57, 59
Active Sequence.....62, 63, 64
Adjust Vial Retainer Arm 77
Agitate Box 92
Agitation Device 90
Agitator Septum.....93
Air Drying 42
Align Syringe Needle 31
Ambient Headspace
Injection..... 5-6, 14, 16, 99, 103, 105, 107
Analysis..... 57
Application Notes 131
AutoControl Module..... 27
AutoLink 60
Automix 65, 66
Automix Routines 65, 66
AutoSampler Timing..... 40, 126

B

Ball Plunger Screws 22
Baseline 57
Build a
Headspace Injection Method 99
SPME Injection Method..... 69, 74
Build a SPME Injection Method 69
Build the Headspace Injection Method..... 99
Build the Liquid Injection Method..... 38

C

Calibration 57
Captive Screws 22
Carrousel..... 5, 7
12 Vial 71
48 Vial 71, 101
Loading 47, 83, 109
Carrousel Label 8, 9, 47, 83

Carrousel Type
12 Vial.....39
48 Vial.....39
Carrousels.....8
Configuration Window . 26-27, 39, 45, 51, 69-70, 101

D

Default Settings.....58
Delete SampleList Entries59
Desorb Time.....71, 86
Diffusion Rate.....89
Downloading a Method.....60

E

EndTime.....40
Equilibrium89

F

Features134
Fiber Holder 16, 69, 72
Fuse.....7, 17, 23

G

Gas Inlet.....32, 103

H

Hardware Configuration27, 38, 69
Headspace Sampling85, 131
Hot Needle Time 16, 44, 102

I

Injection Methods6
Injection Rate 44, 102
Injection Tower Module7, 11, 119
Injector Insert 74, 133
Injector Temperature 85, 127, 136
Installation.....17
Installation Check.....29
Instrument Parameters.....54

L

Large Volume Injection (LVI)111
Leaks48, 121
Liquid Injection 6, 12, 16

- Liquid Injection Method 38
Liquid Injection Operation 47
Liquid Injection Setup 31
Liquid Sampling 16, 28, 85, 101, 133
Lower Air Gap 41, 44, 49
- M**
- Maintenance 117
MANUALn.SEQ 55
MANUALn.SMP 55
Max Errors 54
Method Editor 45, 69, 72, 100
Method File 38, 45, 69
Microvials 15, 16, 43, 50, 123
Mounting Hardware 11, 17
Mounting Plates 17, 19, 22, 75, 106
- N**
- Neat 5, 6, 16, 40, 42, 123, 126
Needle Residence Time 42, 45
New Method Section 38, 70, 100
- P**
- Parts and Supplies 129
Pause Time 40, 44, 102
Performing Injections 51
Pneumatics Module 5, 7, 9, 10, 11, 33, 119
Position Sensor Board 7
Power Switch 7, 9, 20, 23, 26, 36, 51, 74, 106, 118
Pre-drilled Septum 92
Prep Ahead 39, 71, 101
Pressure Regulator 9, 103
- S**
- Sample Vials 50, 83, 84, 109, 123
SampleList 55, 56, 57, 59, 60, 61, 62, 64, 66
Sampling Modes 40, 70, 101
Saturn 28
Semi-volatile Components 89
Septum 13, 37, 43, 47, 49-50, 73, 76, 83, 104,
..... 106, 109, 123, 125-126, 129, 133, 135
Sequence 55, 62, 64
Solvent Flush Sampling 40, 42, 44, 101
Solvent Plug Size 43, 49, 122
Solvent Reservoir Bracket 31, 32, 46
Solvent Reservoirs 9, 47, 104
Solvents 5, 6, 40, 47, 119, 122, 124, 126, 131
Specifications 15, 16
SPME Fibers 130
SPME Holder 13, 73, 130
SPME III Agitation 89
SPME Injection 6, 13, 16, 17, 69, 74
SPME Injection Method 69, 71, 74, 78, 83
SPME Injection Operation 83
Standard 16, 40
Standard Accessory Kit 14
Standard Mode 15, 48, 122
Standard Solvent Flush 6
Star Chromatography Workstation 5, 55, 62, 69
STARTUPn.MTH 64
Status Cable 20
Storage Module 5, 7, 11, 15, 18, 22, 76, 106
Storage Module Latch Release 36, 74, 106
Syringe Carriage 11, 14, 19, 34, 37, 49,
..... 69, 72, 74, 78, 105-106
 Drive Belt, Guide Rods 118
Syringe Carriage Bracket 35
Syringe Carriage Stop Solenoid 36
Syringe/Needle Assembly 12, 31, 34, 46, 105, 120
System Control Configuration 26-27, 51, 56, 63
System Description 7
- T**
- Test Sample 85, 129
Troubleshooting 117
Turntable 7, 8, 9, 83
- U**
- Upper Air Gap 40, 44, 49
Uptake Rate 89
Uptake Speed 16, 40, 42, 44, 102
User-Defined Mode 6, 40
- V**
- Venting 48
Verification 57
Vial Needle Depth 43, 101, 125
Vial Retainer Alignment 25
Vial Retainer Arm 14, 69, 78, 126
 Adjust 77
Vial Sensors 7
Vial Septa 50, 84, 109
Viscous 5, 6, 16, 40, 41, 44, 122, 131
Volatile 5, 6, 16, 40, 123, 131
Volatile Solvents 48
Voltage Selector Card 23
- W**
- Wash Cycle 15, 42, 50, 122, 125
Waste Arm 14, 41, 49, 123

